





17LA3

DE	Gebrauchsanweisung	4
EN	Instructions for use	11
FR	Instructions d'utilisation	18
IT	Istruzioni per l'uso	26
ES	Instrucciones de uso	33
PT	Manual de utilização	41
NL	Gebruiksaanwijzing	48
SV	Bruksanvisning	56
DA	Brugsanvisning	62
NO	Bruksanvisning	69
FI	Käyttöohje	76
CS	Návod k použití	83
SK	Návod na používanie	90
RU	Руководство по применению	97
ZH	使用说明书	105

INFORMATION

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INFORMATION

Last update: 2015-12-10

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide you with important information on the processing of the 17LA3 ankle joints.

2 Product description

2.1 Available sizes

The unilateral system ankle joint is available in four sizes. It can be used either only medially or laterally of the ankle or medially and laterally of the ankle. The choice of system size is based on the patient's body weight and conditions of use.

	Maximum patient weight [kg]			
	17LA3=12/ -T	17LA3=14/ -T	17LA3=16/ -T	17LA3=20/ -T
Unilateral use (lateral or medial)	20	50	85	110
Bilateral use (lateral and medial)	40	80	120	160

INFORMATION

- ▶ When using the product unilaterally, and in case of flexion contracture in the knee or hip > 10° or distinct torsion or valgus/varus instabilities or increased bodily activity, the next higher size must be used.
- ▶ For tubercle seat or valgus/varus malpositions > 10° the product must be fitted bilaterally.

2.2 Components/design (see Fig. 1)

Item	Component	Quantity	Item	Component	Quantity
1	Set screw	2	8	Bushing with collar	2
2	Stop pin (small)	2	9	Screw	1
3	Compression spring	2	10	Flat head screw (TORX®)	1
	Compression spring set (only 17LA3=16/=20)	2			
4	Stop pin	2	11	Axial washer	2
5	Bearing ball	2	12	Flat head screw (Socket head)	2
6	Basic piece	1	13	Stop dummy	2
7	Nut	1			

The scope of delivery does NOT include:

Item	Component	Article number
14	Joint bar	17LS3=12/-T, =14/-T, =16/-T, =20/-T
15	Foot stirrup	17LF3=12, =14, =16, =20
16	Lamination dummy	17LD1=12, =14, =16, =20
17	Shoulder screw (included in scope of delivery of 17LD1=*)	

3 Application

3.1 Indications for use

The system ankle joints are **exclusively** intended to be used for the patient's orthotic fitting of the lower limbs.

In case of unilateral use, the manufacturer recommends that the system ankle joint be processed by means of lamination resin technique or carbon prepreg technology. When used bilaterally, there is no limitation on the processing technology used.

Use of the 17LA3=16-T or 17LA3=20-T system ankle joints in the C-Brace orthotronic mobility system is permitted exclusively with the 17AO1 sensor screw.

3.2 Indications

Partial or total paralysis of the leg muscles.

Indications must be determined by the physician.

4 Safety

4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

4.2 Safety instructions

CAUTION

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service etc.).

CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

CAUTION

Incorrect alignment, assembly or adjustment

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at a temperature above 300 °C (570 °F).

NOTICE

Thermal overloading of product components

Loss of function due to improper thermal treatment

- ▶ Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

5 Preparation for use

5.1 Processing

INFORMATION

The instructions for manufacturing an orthosis with unilateral joint systems are included in the 646T5=9.1 Technical Information. The laminate structure described here is a recommendation of Ottobock and can be adjusted to the needs of the patient.

- > **Required tools and materials:** 646G15 woven carbon stockinette, 17Y106 PVC profile material, 616G12 carbon fibre cloth, 99B81=* PVA bag, 623T3=* Perlon stockinette, 636K8 Plastaband, 616G6 Dacron felt, 633W8 wax, 17LF3=* stirrup, 17LS3=* joint bar, 617H55 C-Orthocryl lamination resin, 743R6 adjustment set for orthotic joints.
- > **Optional:** 17LD1=* lamination dummy, stop dummies
 - 1) Place cut pads or shim layers on plaster model.
 - 2) Pull a layer of perlon stockinette over the plaster model.

- 3) Align the bar parts on the plaster model.
- 4) **When using a lamination dummy:** Fix joint bar and stirrup to alignment adapter (743Y56=4) using the lamination dummy and the enclosed shoulder screw (see fig. 2 and see fig. 4).
- 5) **When not using a lamination dummy:** Fix joint bar and stirrup to alignment adapter (743Y56=4) using the basic piece (no illustration). Place the stop dummies in the spring passages of the basic piece.
- 6) **NOTICE! Make sure that the marking of the foot stirrup is aligned in the neutral position vis-à-vis the marking of the counterpart (see Fig. 3, items 1 and 2).**
- 7) Pull a PVA bag over the plaster model to isolate it.
- 8) Pull a layer of perlon stockinette over the plaster model.
- 9) Fix the finished reinforced bar parts to the plaster model. Carry out the reinforcement of the model in accordance with patient's needs (see fig. 5).
- 10) Seal the lamination dummy or basic piece of the joint with Plastaband and fill hollows with Dacron felt (see fig. 6).
- 11) Pull another layer of perlon stockinette over the plaster model.
- 12) Affix further reinforcement layers for the support surfaces and the foot area, as described in section 9 below.
- 13) Complete reinforcement with a layer of Perlon stockinette and designer fabric.
- 14) Pull on a PVA bag and vacuum.
- 15) Then laminate with C-Orthocryl lamination resin and let it cure.
- 16) Demould the model and cut loose the lamination dummy (16) or basic piece (6) of the system ankle joint.
- 17) **If using without a lamination dummy:** Remove stop dummies from the spring passages.

5.2 Assembly

	Torque [Nm]			
	17LA3=12	17LA3=14	17LA3=16	17LA3=20
Flat head screw (item 10)	4.0	8.0	8.0	12.0

- 1) Insert foot component with the stirrup (15) and axial washers (11) in the basic piece. Use one axial washer per side.
- 2) Make sure that the axial washers prevent direct contact between the metal basic piece and the foot stirrup.
- 3) Lubricate the screw (9) and nut (7) with the 633F7 special lubricant.

- 4) Insert the screw (9) and nut (7) and fasten using the flat head screw (10). Secure the screw connection with Loctite® 241 and tighten to the specified torque value (see table).

5.3 Adaptation

The system ankle joint is adapted to the needs of the patient through the individual combination of the stop pin (4), bearing ball (5), compression spring (3) and stop pin (small) (2) with the set screws (1).

The following combinations are possible:

Combination		Application
Installation position anterior	Installation position posterior	
Stop pin		e.g. ICP or spina bifida
Spring	Stop pin	e.g. in case of distinct hyperextension in the knee joint
Spring stop	Stop pin	e.g. with distinct hyperextension in the knee joint, but with additionally adjustable damped dorsal stop
Spring		e.g. in case of weakness of the lower leg muscles
Stop pin	Spring	Dorsiflexion assist, e.g. in case of peroneal paralysis without own muscular knee joint protection
Spring stop	Spring	e.g. in case of pareses of the lower leg muscles
Spring	Spring stop	e.g. in case of weakness of the lower leg muscles with tendency to hyperextension of the knee joint
Stop pin	Spring stop	Dorsiflexion assist, e.g. in case of pareses of the lower leg muscles without muscular knee joint protection
Spring stop		e.g. in case of pareses in the lower leg muscles with weakness of the knee-protecting muscles

	Torque [Nm]			
	17LA3=12	17LA3=14	17LA3=16	17LA3=20
Countersunk head screws (item 12)	3.6	3.6	7.1	7.1

- 1) Install the bearing ball (5), stop pins (4), compression spring (3) and stop pins (small) (2) with the set screws (1) in the combination required.
- 2) Secure the set screws (1) with Loctite® 241.

- 3) Secure the joint on the joint bar (14) using the flat head screws (12) and secure the screw connection with Loctite® 241 or tighten to the specified torque value (see table).

6 Cleaning

After contact with water containing salt, chlorine or soap, or when they get dirty, the system joints must be cleaned promptly.

- 1) Rinse the system joints with clean fresh water.
- 2) Dry the system joints with a cloth or allow to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

7 Maintenance

INFORMATION

The product may be exposed to increased loads by the patient.

- Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

- Lubricate the springs (3) and the bearing surfaces of the axial washers (11) using the 633F7 special lubricant.

Spare parts are listed in the section "Components/Design" (see fig. 1).

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH