Otto Bock®

QUALITY FOR LIFE


CE

50S3 Malleo Sprint

DE	Sprunggelenkorthese	4
EN	Ankle Orthosis	7
FR	Orthèse pour de maintien pour l'articulation astragalocalcanéenne	11
IT	Cavigliera	14
ES	Órtesis para articulación astrágaloalcalcánea	18
PT	Ortótise para tornozelo	21
NL	Enkelbandage	25
SV	Fotledsortos	28
DA	Ankelleds-ortose	32
NO	Fotleddkortese	35
PL	Orteza stawu skokowego	38
HU	Bokaortézis	41
CS	Ortéza kotníková	45
HR	Ortoza za gležanj	48
TR	Ayak eklemi kemiği ortezi	51
EL	Κνημοποδικός ορθωτικός νάρθηκας	55
RU	Голеноостопный ортез «MalleoSprint»	58
ZH	踝关节矫形器	62

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Explanation of Symbols

English

 CAUTION	Warnings regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.
 INFORMATION	Additional information on the fitting / use.

INFORMATION

Please read these Instructions for Use thoroughly. Follow the safety instructions.

1 Description

1.1 Intended Use

The 50S3 MalleoSprint Ankle Orthosis is to be used **exclusively** for orthotic fitting of the lower limb.

1.2 Indications and Effects

The 50S9 MalleoSprint ankle orthosis is used for functional treatment of ligament instabilities, post-operatively after ligament ruptures, after ankle distortions, as well as for functional-conservative therapy of lateral collateral ligament ruptures. Indication must be determined by the physician.

The orthosis stabilizes the ankle with its special exterior plastic skeleton and the strap guide. Eversion and inversion of the foot joint are restricted without negatively affecting plantar and dorsal flexion of the foot.

1.3 Contraindications

1.3.1 Absolute Contraindications

None.

1.3.2 Relative Contraindications

Skin diseases/injuries, inflammatory symptoms, prominent scars on body parts to which the orthosis is to be applied that are swollen, reddened and feel excessively warm; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the orthosis is to be applied; sensory and circulatory disorders.

1.4 Safety Instructions

INFORMATION

The patient is to be instructed in the proper use / care of the product.

The initial fitting and application of the product must be carried out by trained, qualified personnel only. The daily duration of use and period of application are dependent on medical indication by the physician.

⚠ CAUTION

Risk of injury as a result of improper use. The product is designed for use on one patient only.

An orthosis / support applied too tightly to the body can cause local pressure and, in some cases, even restrict adjacent vessels or nerves. Do not apply the product too tightly. Consult a physician immediately if you experience unusual changes (such as increase in pain).

Improper changes to the product are not permitted.

NOTICE

Risk of damage due to incorrect environmental conditions. The product is not flame-resistant. It must be kept away from flames or other heat sources.

The product should not come into contact with grease or acidic agents, unguents and lotions as this can reduce the period of use—especially of pads.

2 Handling

2.1 Size Selection

Size selection of the orthosis is based on the ankle joint circumference.

2.2 Adapting and Applying the Product

The initial fitting and application of the orthosis must be carried out by qualified personnel.

The orthosis is delivered ready for trial fitting.

1. Open all straps and loosen laces.
2. Pull the orthosis over the foot (fig. 1) and lace up like a shoe (fig. 2). You may also use hooks for quick lacing (fig. 3).
3. Next, apply the stabilizing straps as follows: Lay the outer strap over the instep and pass it from inner (medial) side under the foot (fig. 4) and attach to the corresponding hook and loop closure on the ankle.
4. Next, lay the inner strap over the instep, pass it from the outer (lateral) side under the foot and attach to the corresponding hook and loop closure on the ankle (fig. 5).
5. If full weight-bearing of the foot has been allowed by the physician, the stabilization straps may be re-tightened (fig. 6).

2.3 Instructions for Use and Care

The orthosis must be used in strict conformity with the indication. The physician's instructions must be strictly observed.

Material: Side elements: polyethylene; textile pads: velour/PU; inner side of the textile: PU; hook and loop closures: nylon straps

Cleaning: Use a slightly moistened cloth to clean the outside of the orthosis. Hand-wash the textile elements in 40 °C warm water with a standard mild detergent. We recommend closing the hook and loop closures before washing to keep them functional for a long time. Rinse thoroughly so that all residual detergent is removed. Let air dry. Do not expose to the sun and do not use oven or radiator heat for drying.

If need be, the shoelaces can be replaced.

Disposal after usage: The material used in the product is recyclable and should be disposed of accordingly.

3 Further Usage Restrictions

The product is designed for use on **one** patient only. Re-use of the product is not permitted.

The daily duration of use and period of application are dependent on the medical indication.

4 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the¹ purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

5 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH