3R60 Vacuum

DE Gebrauchsanweisung	5
EN Instructions for Use	15
FR Instructions d'utilisation	26
IT Istruzioni per l'uso	38
ES Instrucciones de uso	49
PT Manual de utilização	60
NL Gebruiksaanwijzing	71
SV Bruksanvisning	81
DA Brugsanvisning	92
NO Bruksanvisning	102
FI Käyttöohje	112
TR Kullanma talimatı	123
RU Руководство по применению	133
ZH 使用说明书	145
KO 사용 설명서	154

1

Aufbaulinie

Alignment
reference line

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INFORMATION

Last update: 2014-06-13

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

1.1 Construction and Function

The 3R60 Vacuum (3R60=VC) is a polycentric prosthetic knee joint with EBS (ergonomically balanced stride) unit and an integrated, mechanical vacuum pump.

The EBS unit permits controlled stance phase flexion up to 15° while simultaneously improving knee stability. This is accomplished by shifting the virtual pivot point (ICR – instantaneous centre of rotation) in the proximal and posterior directions. Stance phase flexion takes place against the resistance of an elastomer spring. Stance phase extension is dampened by hydraulics.

Flexion and extension during the swing phase is controlled by two continuously variable hydraulic valves.

An integrated extension assist supports the swing phase.

The vacuum pump in the prosthetic knee joint is part of a vacuum system that improves the connection between the prosthesis and residual limb. An ischial containment prosthetic socket and other components are needed as well.

1.2 Combination Possibilities

For combination possibilities, please see the catalogue 646K2* or consult the manufacturer.

INFORMATION

Ottobock recommends using the product only in combination with socket systems that distribute the vacuum over the largest possible area. The seal of the prosthetic socket should begin as close underneath the socket entry section as possible. This ensures optimum control of residual limb volume fluctuations.

Ottobock recommends the following socket system:

- ▶ ProSeal Ring (452A1=*) with a textile-free liner (e.g. 6Y81 Skeo Liner) and the 2R119 Vacuum Connection

2 Application

2.1 Area of Application

Area of application according to the MOBIS mobility system:

Recommended for mobility grades **2 and 3** (restricted outdoor walker and unrestricted outdoor walker). Approved for a body weight of up to **125 kg**.

2.2 Environmental Conditions

Allowable environmental conditions

Temperature range for use: -10 °C to +60°C

Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions

Mechanical vibrations or impacts

Perspiration, urine, fresh water, salt water, acids

Dust, sand, highly hygroscopic particles (e.g. talcum)

2.3 Service Life

This prosthetic component has been tested by the manufacturer for 3 million load cycles according to ISO 10328. Depending on the patient's activity level, this corresponds to a service life of 3 to 5 years.

2.4 Qualification

The product may only be integrated in the prosthesis by orthopaedic technicians who have completed training for the Harmony system and obtained corresponding certification.

3 Safety

3.1 Explanation of Warning Symbols

⚠ WARNING Warning regarding possible serious risks of accident or injury.

⚠ CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

3.2 General Safety Instructions

⚠ WARNING

Operation of motor vehicles

Risk of accidents due to restricted body functions

- ▶ Observe the applicable legal and insurance regulations for the operation of motor vehicles and have your driving ability examined and certified by an authorised agency.
- ▶ **Inform the patient.**

⚠ CAUTION

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

⚠ CAUTION

Excessive strain on the product

Fall due to breakage of load-bearing components

- ▶ Use the prosthesis components according to the MOBIS classification (see section "Area of Application").
- ▶ **Inform the patient.**

⚠ CAUTION

Unallowable combination of prosthesis components

Fall due to breakage or deformation of the product

- ▶ Only combine the product with prosthesis components that are approved according to the section "Combination Possibilities".
- ▶ Based on the instructions for use of the prosthesis components, verify that they may be combined with each other.

CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if applicable (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).
- ▶ **Inform the patient.**

CAUTION

Exceeding the service life and reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded (see section "Service Life").
- ▶ Only use the product for a single patient.
- ▶ **Inform the patient.**

CAUTION

Reaching into the area of the joint mechanism

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.
- ▶ **Inform the patient.**

CAUTION

Use of talcum

Falling, damage to the prosthetic knee joint due to lack of lubrication

- ▶ Do not use talcum on the prosthetic knee joint or other prosthesis components.
- ▶ **Inform the patient.**

CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check a damaged product for proper functionality and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take suitable measures if required (e. g. repair, replacement, inspection by the manufacturer's customer service etc.).
- ▶ **Inform the patient.**

Signs of changes in or loss of functionality during use

Noticeable functional changes may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc.

4 Scope of Delivery

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering according to their respective symbols as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

3R60 Vacuum					
Fig.		Item no.	Quantity	Designation	Article number
–	■	–	1	Instructions for use	647G873
3	–	–	1	3R60 Vacuum	3R60=VC
3	■	–	2	Tube	4X145
–	■	–	1	Filter	4X253
3	■	–	1	Cosmetic Exhaust Flange	4Y383
–	■	–	1	Adjustment Wrench	710H10=2x3

5 Preparation for Use

 **CAUTION**

Incorrect alignment, assembly or adjustment
Injuries may result from incorrectly assembled, incorrectly adjusted or damaged components

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.
- ▶ **Inform the patient.**

5.1 Preparing the prosthetic socket

 **CAUTION**

Using a wrong prosthetic socket
Contusions and pressure sores due to wrong socket fit

- ▶ Use the vacuum pump only in combination with an ischial containment prosthetic socket with corresponding end load.

Mounting the Vacuum Connector

- ▶ Mount the 2R119 Vacuum Connector on the prosthetic socket according to the applicable instructions for use.

5.2 Bench Alignment

Bench alignment process	
+ = anterior placement / – = posterior placement (in relation to the alignment reference line)	
Item	see Fig. 1, see Fig. 2
	Required tools and materials: 662M4 Goniometer, 743S12 Heel Height Measuring Apparatus, 743A80 50:50 Gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly)
	Position prosthetic components in accordance with the following values in the alignment apparatus and install them in accordance with the instructions.

Bench alignment process	
+ = anterior placement / – = posterior placement (in relation to the alignment reference line)	
Item	see Fig. 1, see Fig. 2
①	a–p positioning, middle of the prosthetic foot to the alignment reference line: +30 mm
②	Heel height: effective heel height (x) + 5 mm Exterior foot rotation: approx. 5°
③	Positioning of the alignment reference point (front upper axis of the prosthetic knee joint) in the alignment apparatus: MTP-floor measurement + 20 mm a–p positioning of the alignment reference point to the alignment reference line: 0 mm External rotation of the prosthetic knee joint: approx. 5°
④	Join the prosthetic foot with the prosthetic knee joint using the tube adapter. Observe the adapter instructions for use during adaptation and installation.
⑤	Mark 2 centre points laterally on the prosthetic socket: 1. At the level of the trochanter major 2. In the distal region Join the two points with a line.
⑥	Position the prosthetic socket so its proximal point lies on the alignment reference line. Set the socket flexion to 3° to 5° . Take the individual situation of the patient into account.
⑦	Join the prosthetic socket and the prosthetic knee joint using the adapters you selected. Be sure to follow the adapter instructions for use during adjustment and installation.

5.3 Static Alignment

Static alignment process	
+ = anterior placement / – = posterior placement (in relation to the alignment reference line)	
Item	see Fig. 1
	Required materials and tools: 743L100 L.A.S.A.R. Posture
①	To determine the load line, position the patient on the L.A.S.A.R. Posture as follows: - Prosthesis side on force measurement plate (sufficient load: > 35% body weight) - Other side on height compensation plate
②	Optimise the static alignment solely by changing the plantar flexion. a–p positioning of the front lower axis of the prosthetic knee joint to the load line: -10 mm

5.4 Installing the hoses

The hoses can be installed after completing the bench alignment and static alignment.

5.4.1 Installing the intake hose

The intake hose runs from the suction connector of the vacuum pump to the vacuum connector in the prosthetic socket. The filter inserted in the hose protects the vacuum pump from contamination.

Installing the filter

- 1) Cut the hose at the appropriate place (see Fig. 11).
- 2) Connect the filter to the two hose ends (see Fig. 12).

Installing the intake hose

- 1) Align the intake hose in the proximal direction and connect it to the suction connector of the vacuum pump (see Fig. 4).
- 2) **INFORMATION: If a rotation adapter is used, the intake hose has to be longer.**
Run the intake hose to the vacuum connector of the prosthetic socket and connect it (see Fig. 5).

5.4.2 Installing the outlet hose

CAUTION

Improper installation of the outlet hose

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Run the outlet hose in such a way that the hose opening is not directed towards parts of the body or into the cosmetic cover.

The suctioned air and fluids are drained through the outlet hose. Fluids must not be drained onto metal components of the modular prosthesis. Increasing the length of the outlet hose lowers the noise level of the vacuum pump.

- 1) Align the outlet hose in the distal direction and connect it to the corresponding connector on the vacuum pump (see Fig. 6).
- 2) **Optional:** Clamp the outlet hose in the hose clamp of the prosthetic knee joint (see Fig. 7).
- 3) **If a cosmetic foam cover is used:** Make an opening in the cosmetic foam cover and install the 4Y483 Cosmetic Exhaust Flange.
- 4) **If no cosmetic cover is used:** Attach the outlet hose to the prosthesis in a suitable position.
- 5) Trim any excess outlet hose.

5.5 Dynamic Trial Fitting

CAUTION

Dynamic trial fitting without parallel bars

Fall due to lack of patient experience or functional changes

- ▶ Perform the dynamic trial fitting using parallel bars when the prosthesis is first used and also after any change is made to the prosthesis.

CAUTION

Reaching into the area of the joint mechanism

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.
- ▶ **Inform the patient.**

5.5.1 Testing the vacuum system

The functionality of the vacuum system is verified in the course of dynamic trial fitting. A vacuum gauge is installed in the system for this purpose.

> **Required materials and tools:**

755Z39 Vacuum Gauge (with T-piece)

- 1) Disconnect the hose from the suction connector of the prosthetic knee joint.
- 2) Connect the vacuum gauge to the suction connector using the T-piece.
- 3) Connect the hose to the vacuum gauge.
- 4) Attach the vacuum gauge to the prosthesis so that the patient is able to move freely.
- 5) Conduct the dynamic trial fitting. Check the vacuum during trial fitting.
- 6) Check the system if there are any leaks (see Page 24).
- 7) Remove the vacuum gauge once testing is complete.

5.5.2 Setting flexion and extension damping

The desired gait pattern is achieved by adjusting the damping behaviour of the prosthetic knee joint for flexion and extension in the swing phase.

INFORMATION

The valves for flexion and extension are located on the left and right sides of the prosthetic knee joint.

Valve	Symbol	Factory setting (damping)	
Flexion valve	(F)	Centre position	Medium
Extension valve	(E)	Stop –	Low

> **Required tools and materials:**

710H10=2x3 Adjustment Wrench

- 1) Use the adjustment wrench to set the flexion resistance (see Fig. 8).
- 2) Use the adjustment wrench to set the extension resistance (see Fig. 9).

Situation	Measures	Setting
Lower leg swings through too far	Increase damping	Turn flexion valve to the right (+)
Lower leg does not swing through sufficiently	Reduce damping	Turn flexion valve to the left (–)
Lower leg impact too hard on the extension stop	Increase damping	Turn extension valve to the right (+)
Lower leg does not reach full extension before heel strike	Reduce damping	Turn extension valve to the left (–)

5.5.3 Adjusting stance phase flexion

Stance phase flexion takes place upon heel strike while walking and is controlled by the EBS unit. A nut is used for pre-compression of the elastomer spring to control the response behaviour. Hydraulic damping for stance phase extension cannot be changed. The transition to the flexion movement of the swing phase begins when the ball of the foot makes contact.

INFORMATION

The adjustment nut of the EBS unit is found on the back of the prosthetic knee joint below the elastomer block.

> **Required tools and materials:**

710H10=2x3 Adjustment Wrench

- > The patient has been instructed in the functionality of the EBS unit.
- 1) Check the behaviour of the EBS unit during trial walking.

2) **INFORMATION: Reducing the pretension increases knee safety while increasing the pretension decreases knee safety.**

Adjust the pretension of the EBS unit in order to achieve the desired stance phase flexion behaviour (see Fig. 10). Ensure that the pivoting mount does not permanently touch the end stop while walking.

Situation	Conclusion	Measures
Strong stance phase flexion, patient experiences yielding until the pivoting mount reaches the end stop.	EBS unit too soft (pretension too low)	Turn the adjustment nut to the right (+)
Little stance phase flexion, yielding is not sufficient for the patient.	EBS unit too firm (pretension too high)	Turn the adjustment nut to the left (-)

6 Use

CAUTION

Overheating of the hydraulics due to excessive strain

Burns, injuries caused by falling due to changes in functionality and damage to prosthetic components

- ▶ Do not touch overheated prosthetic components.
- ▶ In case of changes in functionality, stop all activities and allow the overheated prosthetic components to cool down.
- ▶ In case of changes in functionality, reduce strain on the prosthetic component in order to allow cooling.
- ▶ In case of overheating or changes in functionality, have the prosthetic component inspected for damage by a service centre authorised by the manufacturer.
- ▶ **Inform the patient.**

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).
- ▶ **Inform the patient.**

6.1 Using the vacuum pump

The vacuum pump is operated mechanically while walking with the prosthetic knee joint. A vacuum is created in the prosthetic socket during flexion while the medium being conveyed is exhausted through the discharge valve during extension. An elevated vacuum is reached in the sealed socket environment after several pump cycles. A vacuum can be generated even before walking by flexing and extending the prosthetic knee joint while standing.

6.2 Flushing the vacuum pump

CAUTION

Improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Clean the product regularly.

A reduced vacuum indicates contamination. In case of contamination, the system has to be flushed to clean it. This flushes the vacuum pump and valves. The flushing process can also help if the system is plugged.

> **Required materials:**

Vessel with 30 ml to 60 ml distilled water, one collecting vessel

Preparations required:

The prosthesis is not worn for this process

- 1) Pour the distilled water into the prosthetic socket.
- 2) Position the vessel under the outlet hose.
- 3) Flex and extend the prosthetic knee joint to pump the water through the system.
- 4) Continue performing these movements until the full volume of water has been pumped out of the system.
- 5) Clean the interior surface of the socket.

7 Cleaning

CAUTION

Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Only clean the product with the approved cleaning agents.
- ▶ Only disinfect the product with the approved disinfectants.
- ▶ Observe the instructions for cleaning and care.

Inform the patient.


- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

8 Maintenance

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8.1 Checking for leaks


Remove the hose between the vacuum connector of the prosthetic socket and the suction connector of the prosthetic knee joint. Use the T-piece to connect the vacuum gauge to the vacuum connector and suction connector. Cover the entire vacuum connector on the inside of the socket with a small piece of paper and airtight tape (e.g. Coroplast). Flex and extend the prosthetic knee joint until a constant vacuum of at least **300 mbar** is reached. Check whether the vacuum level is maintained.

Problem: Proximal seal

Solution: Check the proximal seal; seal or replace as necessary.


Remove the adhesive tape and piece of paper. Cover the hole in the vacuum connector with a smaller piece of paper and tape. Flex and extend the prosthetic knee joint until a vacuum of at least **300 mbar** is reached. Check whether the vacuum level is maintained.

Problem: Vacuum connector sealing

Solution: Check the vacuum connector sealing for leaks; seal or replace as necessary.


Connect the vacuum gauge to the hose of the prosthetic knee joint. Flex and extend the prosthetic knee joint until a vacuum of at least **300 mbar** is reached. Check whether the vacuum level is maintained.

Problem: Vacuum connector insert

Solution: Replace the vacuum connector insert.

Method of resolution: Replace the hose between the vacuum gauge and prosthetic knee joint. Flex and extend the prosthetic knee joint until a vacuum of at least **300 mbar** is reached. Check whether the vacuum level is maintained.


Problem: Hose

Solution: Replace the hose between the vacuum connector of the prosthetic socket and the suction connector of the prosthetic knee joint.

Problem: Valves or vacuum pump

Solution: Remove the valves of the 3R60 Vacuum; clean or replace them. Connect the vacuum gauge to the suction connector of the prosthetic knee joint. Flex and extend the prosthetic knee joint until a vacuum of at least **300 mbar** is reached. Check whether the vacuum level is maintained.

Problem: Vacuum pump

Solution: Please contact the manufacturer's local service representative and send in the defective prosthetic knee joint.

✓ = The vacuum does not decrease within 20 seconds.

✗ = The vacuum decreases within 20 seconds.

9 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

10 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

11 Technical data

Article number	3R60=VC
Weight [g]	890
Proximal system height up to the alignment reference point [mm]	1
Distal system height [mm]	173
Proximal connection	Pyramid adapter
Distal connection	Pyramid adapter
Max. flexion angle [°]	173
Max. body weight [kg]	125
Mobility grade	2 to 3

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Otto Bock HealthCare GmbH