





Kenevo Protector 4X840

DE Gebrauchsanweisung (Fachpersonal)	4
EN Instructions for use (qualified personnel)	9
FR Instructions d'utilisation (Personnel spécialisé)	14
IT Istruzioni per l'uso (Personale tecnico specializzato)	19
ES Instrucciones de uso (Personal técnico especializado)	24
PT Manual de utilização (Pessoal técnico)	30
NL Gebruiksaanwijzing (Vakmensen)	35
SV Bruksanvisning (Fackpersonal)	40
DA Brugsanvisning (Faguddannet personale)	45
NO Bruksanvisning (Fagpersonell)	50
FI Käyttöohje (Ammattihenkilöstö)	55
PL Instrukcja użytkowania (Personel fachowy)	60
HU Használati utasítás (szakszemélyzet)	65
CS Návod k použití (Odborný personál)	70
SK Návod na používanie (Odborný personál)	75
TR Kullanma talimatı (Uzman personel)	80
RU Руководство по применению (Квалифицированный персонал)	85

INFORMATION

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INFORMATION

Last update: 2015-07-08

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- Please keep this document in a safe place.

The "Kenevo Protective Cover" is referred to simply as the protective cover / product below.


These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product Description

2.1 Design

The product consists of the following components:

1. Protective cover main unit
2. Upper closure (one side always stays on the main unit)
3. Lengthwise closure
4. Flexible cover cap

2.2 Function

The protective cover provides a means of enclosing the Kenevo electronic knee joint, including the tube adapter, to enhance its appearance and function. The housing of the knee joint is protected against external influences by the protective cover, e.g. scratching and surface contamination. The protective cover is not required for the proper and safe functioning of the joint.

2.3 Combination Possibilities

Prosthetic components that are suited for combination with the product are listed below.

Knee Joint Combination Possibilities

Designation	Reference number
Kenevo	3C60/3C60=ST

3 Application

3.1 Indications for use

The housing of the knee joint is protected against external influences by the protector, e. g. scratching and surface contamination.

3.2 Area of Application

The product is intended **exclusively** for wearers of the Kenevo knee joint.

3.3 Conditions of use

The product was developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

3.4 Qualification

The fitting of patients with the product may only be carried out by qualified personnel (e.g. a prosthetist).

4 Safety

4.1 Explanation of Warning Symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 CAUTION
The heading describes the source and/or the type of hazard The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible: > E.g.: Consequence 1 in case of failure to observe the hazard > E.g.: Consequence 2 in case of failure to observe the hazard ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General Safety Instructions

 NOTICE
Improper assembly/disassembly of the product Damage to the product.

- Observe the process steps, instructions and required tools in the section "Preparation for use".
- Instruct the patient in the proper disassembly/ assembly of the product.

NOTICE

Independent user changes or modifications made to the product

Damage to the product.

- Manipulations to the product other than the tasks described in these instructions for use are not permitted.

NOTICE

Mechanical damage to the product

Change in or loss of functionality due to damage.

- Use caution when working with the product.
- If the product is damaged, check it for proper function and readiness for use.
- Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

5 Scope of Delivery

- 1 pc. Kenevo Protective Cover main unit
- 1 pc. 4P840=R lengthwise closure
- 1 pc. 4P860=U upper closure
- 1 pc. 4P8 flexible cover cap
- 1 pc. 647G1139 instructions for use (qualified personnel)
- 1 pc. 647G1138 instructions for use (user)

6 Preparation for Use

Required Tools

Designation	Reference number
Vibrating/band saw	-
Knife	-
Loctite® 241 Thread Locking Compound	636K13
Crepe Adhesive Tape	627B6=15
Torque Wrench	-
Tape Measure	-

6.1 Shortening protective cover

- 1) Remove the prosthesis from the residual limb.
- 2) Check whether the closures on the protective cover are properly fastened and that the protective cover is not twisted.
- 3) Lay the protective cover next to the prosthesis.
- 4) Measure the length from the centre of the knee to the lower edge of the tube adapter on the prosthesis (see fig. 1).
- 5) Mark the measured length at at least three places all around the protective cover (see fig. 2).
- 6) **INFORMATION: For tall patients, shortening the protective cover may not be required.**
- 7) Apply crepe adhesive tape proximally all around at the mark.

- 8) Use a vibrating or band saw to shorten the protective cover along the bottom edge of the crepe adhesive tape (see fig. 3).

INFORMATION: Make sure that the closures are not shifted while sawing.

- 9) Deburr the cut edge.
- 10) Using suitable tools, remove the profile rib on the interior of the protective cover (see fig. 4).

NOTICE

Collision between protective cover and prosthetic foot

Damage to or functional impairment of the protective cover or prosthetic foot

- ▶ Carry out a function test with the user wearing the protective cover.
- ▶ If the protective cover collides with the prosthetic foot, shorten the protective cover again.

6.2 Protector Installation

- 1) Remove the prosthetic foot from the tube adapter.
- 2) Slide the flexible cover cap over the tube adapter (see fig. 5).
- 3) Install the prosthetic foot on the tube adapter (see fig. 6).
- 4) Apply Loctite® 241 to the two set screws and screw them in.

⚠ CAUTION

Incorrectly secured screws

Falling due to breakage of load-bearing components caused by screw connections coming loose.

- ▶ The set screws in the tube adapter must be secured before they are tightened to the specified torque.
- ▶ See the instructions for use for the knee joint for the specified torque of the set screws.

- 5) Remove the crepe adhesive tape from the protective cover after final shortening.
- 6) Open or remove the closures.
- 7) Spread open the protective cover (see fig. 7).
- 8) Apply the protective cover to the knee.

INFORMATION: The retainers of the protective cover have to engage on the frame element (see fig. 8).

- 9) Fasten the upper closure (see fig. 9) and lengthwise closure (see fig. 10) on the protective cover.
- 10) Fit the flexible cover cap to the distal end of the protective cover main unit. Align the lengthwise grooves in the extension of the lengthwise closure on the dorsal side of the protective cover (see fig. 11).

NOTICE

Loss of stability due to inadequately installed retainers

No hold of the protective cover on the knee joint.

- ▶ Ensure that all retainers and closures are properly engaged and/or fastened.

6.3 Functional Test

After assembly, the following points have to be checked:

- Has any collision between the protective cover and prosthetic foot been eliminated?
- Is the function of the rotation adapter impaired?

Material may be removed from the protective cover to improve the functionality.

INFORMATION: If material is removed, the stability/flexibility of the protective cover on the knee joint may be impaired.

7 Maintenance

- During the standard inspection of the prosthesis components, also inspect the product.

7.1 Cleaning and Care

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- Only clean the product with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

8 Legal Information

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as class I accessories according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

8.3 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

9 Appendices

9.1 Symbols Used

9.1.1 Symbols on the product

Legal manufacturer

Declaration of conformity according to the applicable European directives

LOT PPPP YYYY WW

Lot number

9.2 Technical data

Ambient conditions	
Storage and transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage and transport without packaging	-25 °C/-13 °F to +70 °C/+158 °F max. 93% relative humidity, non-condensing
Operation	-10 °C/+14 °F to +40 °C/+104 °F max. 93% relative humidity, non-condensing

General	
Reference number	4X840
Product service life	Wear part subject to normal wear and tear
Protective cover weight (with closures)	380 g
Flexible cover cap weight	11 g

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Otto Bock Healthcare Products GmbH