






28L10


DE	Gebrauchsanweisung (Fachpersonal)	5
EN	Instructions for Use (Qualified Personnel)	9
FR	Instructions d'utilisation (Personnel spécialisé)	13
IT	Istruzioni per l'uso (Personale tecnico specializzato)	17
ES	Instrucciones de uso (Personal técnico especializado)	22
PT	Manual de utilização (Pessoal técnico)	26
NL	Gebruiksaanwijzing (Vakmensen)	30
SV	Bruksanvisning (Fackpersonal)	34
DA	Brugsanvisning (Faguddannet personale)	38
NO	Bruksanvisning (Fagpersonell)	42
PL	Instrukcja użytkowania (Personel fachowy)	46
HU	Használati utasítás (szakszemélyzet)	51
CS	Návod k použití (Odborný personál)	55
SL	Navodila za uporabo (Strokovno osebje)	59
BG	Инструкция за употреба (Специалисти)	63
TR	Kullanma talimatı (Uzman personel)	67
EL	Οδηγίες χρήσης (Τεχνικό προσωπικό)	71
RU	Руководство по применению (Квалифицированный персонал)	76
JA	取扱説明書 (有資格担当者)	81
ZH	使用说明书 (专业人员)	85


6


7


8


9


10


11


12


13

14

15

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1 Foreword

INFORMATION

- Date of the last update: 2013-12-16
- Please read this document carefully.
 - Follow the safety instructions.

These instructions for use provide you with important information for properly applying and removing the **Tübingen Hip Flexion and Abduction Orthosis according to Prof. Bernau** and for correctly adapting it to the patient.

2 Product Description

2.1 Components (see Fig. 1)

Item	Designation	Item	Designation
1	Shoulder harness with terry cloth covers	6	Push-pull device for spreader bar (red)
2	Hook-and-loop closure with caterpillar graphic	7	Spreader bar with ratchet mechanism
3	Front closures (white)	8	Back closures (red)
4	Beaded cords (green/white)	9	Lower closures (red)
5	Thigh shells	–	2 terry cloth covers and 1 hook-and-loop closure are included as replacements

2.2 Design

The thigh shells (5) and the shoulder harness (1) are joined together by green/white beaded cords, which are held by the white closures (3) and the red closures (8, 9). The hook-and-loop closure (2) and the white closures (3) on the front of the shoulder harness (1) are used when putting on and taking off the Tübingen Hip Flexion Orthosis.

The attending physician can change the length of the beaded cords (4) with the red closures (8, 9) on the shoulder harness (1) and on the thigh shells (5). This establishes the hip flexion position and can be adjusted as the infant grows.

After loosening the push-pull device (6) on the spreader bar (7), the two sections of the spreader bar (7) can be simultaneously adjusted to achieve the desired hip abduction. The numbered scale on the bottom of the spreader bar (7) makes it easy to duplicate the setting.

3 Application

3.1 Intended Use

The 28L10 Tübingen Hip Flexion Orthosis is intended **exclusively** for the treatment of hip dysplasia in infants (< 1 year).

3.2 Indications

Treatment of hip dysplasia in infants (< 1 years): type IIa, IIb, IIc hips, stable according to Graf.

Indications must be determined by the physician.

3.3 Effects

Flexion of the hip joints in excess of 90° with controlled medium abduction is the best prerequisite for fast maturation of hip joints suffering from delayed development. The Tübingen Hip Flexion Orthosis has proven to be highly effective for establishing this hip position. The thigh shells (5) are connected to the shoulder harness (1) with 2 beaded cords (4) used to establish the reproducible hip flexion setting. The desired degree of abduction is adjustable using a spreader bar (7).

3.4 Contraindications

Consult the attending physician if any of the following occur: skin diseases/injuries, inflammation, reddening or warming in the area where the orthosis is fitted.

3.5 Additional Restrictions for Use

The Tübingen Hip Flexion Orthosis was designed for use on **one** infant (< 1 year). Re-use of the product is not permitted.

4 Safety

4.1 Explanation of Warning Symbols

Warning regarding possible risks of accident or injury

Warning regarding possible technical damage

4.2 General Safety Instructions

Incorrect application of the shoulder harness (1)

Shortness of breath, constriction of the respiratory system

- ▶ Apply the shoulder harness (1) correctly (see the section "Application"). The shoulder harness (1) must be applied so that the hook-and-loop closure (2) with the caterpillar is positioned on the infant's chest.
- ▶ Instruct the parents in correct application (see information for parents).

Incorrect application of the orthosis

Malposition of the legs, impairment of hip development

- ▶ Apply the orthosis correctly (see the section "Application").
- ▶ **Instruct the parents in correct application (see information for parents).**
- ▶ **Inform the parents that changes to the red closures (8, 9) and spreader bar (7) are not permitted.**

5 Usage

5.1 Size Selection

The Tübingen Hip Flexion Orthosis is available in 3 sizes. It is selected according to the age of the infant (see sizing table).

5.2 Initial Fitting (by the physician)

The custom fitting of the orthosis by the attending physician is a key factor for good acceptance and successful treatment. The orthosis setting is symmetrical upon delivery. Individual size adjustment is only possible while the orthosis is being worn.

INFORMATION

While putting on the orthosis for the first time, demonstrate the individual application steps for the parents using the included information for parents. The parents must only use the hook-and-loop closure (2) and the white clamp closures (3) to put on and take off the orthosis.

Application

- 1) Open the hook-and-loop closure (2) and the white closures (3), and take the beaded cords (4) out of the closures (3). Spread out the orthosis on the changing table. Position the thigh shells (5) with the beaded cords (4) to the sides of the child (see Fig. 2).
- 2) Apply the shoulder harness (1) around the neck from behind and fasten it loosely on the chest with the hook-and-loop closure (2).

INFORMATION: Ensure that space of at least 2 finger widths is maintained between the hook-and-loop closure (2) and the neck.

- 3) Grasp the legs by the thighs and move them to hip flexion of more than 90° with abduction. To maintain hip flexion, support the child's feet on your abdomen or upper body. Place the thighs into the thigh shells (see Fig. 3).
- 4) Insert the beaded cords (4) into the white closures (3) (see Fig. 4).
- 5) Fasten the white closures (3) (see Fig. 5).

Adjusting the Flexion Position of the Hip Joint

The distance between the shoulder harness (1) and thigh shells (5) determines hip flexion and can be adjusted using the front white closures (3) and the back red closures (8). **Ensure that the hip joints are flexed in excess of 90°** (see Fig. 6) and that the tension of the beaded cords (4) is even.

- 1) If flexion in the hip joint is too low, symmetrically shorten the beaded cords (4) on the front white closures (3). The green-white colour of the beads makes it easier to obtain a symmetrical setting.
- 2) If an additional adjustment on the back red closures (8) is required, consult the section "Adjusting the Orthosis for Growth".

If the tension on the front and back sections of the beaded cords is different or the front hook-and-loop closure (2) is too close to the neck, change the position as follows:

- 1) Open the lower closure (9) (see Fig. 14).
- 2) Release the beaded cord (4) from the lower closure (9): hold the beaded cord with the middle and index fingers, and pull to release the secured bead.
- 3) Slide the thigh shells along the beaded cord (4) until the distance between the shoulder harness (1) and the neck is 2 finger widths on all sides and the thigh shells (5) are evenly applied. The spreader bar (7) must lie symmetrically between the thighs.
- 4) Secure the beaded cord (4) in the lower closure (9): push the beaded cord into the bead holder until the secured bead engages with a click.
- 5) Fasten the lower closure (9).

Documenting Settings

- For reproducible documentation of the settings for the parents, cut off the beaded cords 3 beads above the front clamp closure (3) (see Fig. 7) and/or mark the bead being inserted in the clamp closure with a waterproof marker.

Abduction Setting

- 1) Slide the push-pull device (6) on the spreader bar (7) all the way to the left (see Fig. 8).
- 2) Slightly tip up the spreader bar (7) through slight interior rotation of the infant's legs. Extend the spreader bar (7) no further than the marking (see Fig. 9) and engage it at the chosen abduction setting (see Fig. 10).
- 3) To lock, slide the push-pull device (6) all the way to the right against the stop until the push-pull device engages with audible feedback (see Fig. 11).

5.3 Adjusting the Orthosis for Growth

Depending on the course of treatment and the growth of the child, the Tübingen Hip Flexion Orthosis must be checked and adjusted by the physician. In order to do so, lengthen the beaded cords on the back closures (8). The extra length of the beaded cords is found behind the red closures (8) on the shoulder harness (1). The green-white colour of the beads helps obtain a symmetrical setting for the beaded cords.

- 1) Release the beaded cord (4) from the back closure (8): hold the beaded cord with the middle and index fingers, and pull hard until the secured bead is released. This is meant to be tight since it prevents the beaded cord from coming loose accidentally.
- 2) Lengthen the beaded cord (4).
- 3) Secure the beaded cord (4) in the back closure (8): hold the beaded cord with the middle and index fingers, and pull so that the secured bead engages with a click (see Fig. 12). Insert the excess beaded cord (4) in the guide groove at a 90° angle and secure it in the holder under the terry cloth cover (see Fig. 13).

Proceed the same way on the opposite side and then check to ensure the setting is symmetrical.

5.4 Cleaning

The orthosis is corrosion-protected and washable. The terry cloth covers for the shoulder harness (1) and the hook-and-loop closure (2) may be removed and washed. Use the supplied replacement covers and hook-and-loop closure during cleaning.

- 1) Press on the sides of the closure with your thumbs, unhook the hook-and-loop closure (2) from the plastic loops (see Fig. 15) and remove the hook-and-loop closure.
- 2) Pull the terry cloth covers off the shoulder harness (1).
- 3) Wash the terry cloth cover and the fastened hook-and-loop closure (2) in warm water at **40 °C/104 °F** using standard mild detergent. **Recommendation:** wash in a mesh laundry bag. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).
- 5) Slide the terry cloth covers onto the straps of the shoulder harness (1) and engage the hook-and-loop closure (2) in the plastic loops. Ensure that the hook-and-loop closure engages properly.

5.5 Removal

INFORMATION

To apply and remove the orthosis, the parents must only use the hook-and-loop closure with the caterpillar (2) and the white closures (3).

Removal

- 1) Open the hook-and-loop closure (2) and the white closures (3).
- 2) Take out the beaded cords (4) and lift the child out of the orthosis.

6 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

6.3 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

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Otto Bock HealthCare GmbH