
50P13 Manu Sensa

DE Gebrauchsanweisung	3
EN Instructions for Use	7
FR Instructions d'utilisation	11
IT Istruzioni per l'uso	16
ES Instrucciones de uso	21
PT Manual de utilização	25
NL Gebruiksaanwijzing	30
SV Bruksanvisning	34
DA Brugsanvisning	39
NO Bruksanvisning	43
FI Käyttöohje	47
PL Instrukcja użytkowania	52
HU Használati utasítás	56
CS Návod k použití	61
HR Upute za uporabu	65
SK Návod na používanie	70
TR Kullanma talimatı	74
EL Οδηγίες χρήσης	78
RU Руководство по применению	83
ZH 使用说明书	88

	Seite / Side	Größe / Size	Handgelenkumfang / Wrist circumference	
	L/R	XS	14 – 15 cm	5.5 – 5.9 inch
	L/R	S	15 – 16 cm	5.9 – 6.3 inch
	L/R	M	16 – 17 cm	6.3 – 6.7 inch
	L/R	L	17 – 18 cm	6.7 – 7.1 inch
	L/R	XL	18 – 19 cm	7.1 – 7.5 inch
	L/R	XXL	19 – 20 cm	7.5 – 7.9 inch

Material	Gestrick / Knit Fabric: PA, EL
	Frotteebezüge / Terry Cloth Covers: Baumwolle / Cotton, PES
	Kunststoffschiene / Plastic Splint: PA
	Handgelenkgurt / Wrist Strap: PA, EL
	Pelotte / Pad: viskoelastisches Technogel® / viscoelastic Technogel®

1 Vorwort

INFORMATION

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INFORMATION

Date of the last update: 2012-08-15

- Please read this document carefully.
- Follow the safety instructions.

These instructions for use provide important information on the fitting and application of the 50P13 Manu Sensa wrist support.

2 Application

2.1 Intended Use

The support is intended **exclusively** for orthotic fittings of the wrist joint and **exclusively** for contact with intact skin.

The support must be used according to the indications and following the instructions of a physician.

2.2 Indications

- Irritation (e. g. tendovaginitis, tendomyopathy, minor joint effusion, post-traumatic [distorsion], post-operative, reactive irritation and arthralgia)
- Rheumatoid arthritis
- Carpal tunnel syndrome
- Degenerative wrist disorders (osteoarthritis)
- Feelings of instability (e.g. snapping ulna)

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases / injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support / orthosis will be applied; sensory and circulatory disorders in the body region being treated and in the area of the fingers.

2.4 Effects

The form-knitted support provides targeted wrist compression and stabilisation. The support promotes sensorimotor functions and activates muscular stabilisation. The anatomically pre-shaped and individually adaptable volar plastic splint keeps the hand in an ideal position and limits wrist movement. The pad in the radial section of the wrist joint enhances the micro-massage and promotes the resorption of oedemas, haematomas and effusions.

3 Safety

3.1 Explanation of Warning Symbols

 CAUTION	Warning regarding possible risks of accident or injury
 NOTICE	Warning regarding possible technical damage

3.2 General Safety Instructions

 CAUTION
Using the product without following the instructions for use Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

CAUTION

Contact with heat, embers or fire

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.
- ▶ **Inform the patient.**

CAUTION

Reuse on other patients and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

NOTICE

Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.
- ▶ **Inform the patient.**

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Size Selection

- 1) Determine the wrist circumference.
- 2) Select the size of the support (see sizing table).

4.2 Adaptation and Application

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the support is applied properly and fits correctly.
- ▶ **Inform the patient.**

NOTICE

Improper moulding or application

Damage to the product due to overloading of the material and improper fit of the product due to breakage of load-bearing components

- ▶ The product may only be moulded by qualified personnel.
- ▶ Do not make any improper changes to the product.
- ▶ Always apply the product according to the information in the instructions.
- ▶ **Inform the patient.**

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear and damage prior to each use.
- ▶ Advise the patient that the product must not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

The volar plastic splint is anatomically pre-shaped. In the course of the fitting, the plastic splint can be heated and then individually shaped by the orthotist. The support is delivered closed and ready for trial fitting.

- 1) Open the elastic strap.
- 2) Insert the plastic splint into the pocket on the inside of the support (see Fig. 1). Ensure that the lower end of the plastic splint is in the pocket intended for this purpose.

- 3) Pull the support over the hand, taking into account the opening provided for the thumb (see Fig. 2).
- 4) The elastic strap can be used for additional stabilisation if required (see Fig. 4).
- 5) When providing the support to the patient, the proper fit and positioning of the support must be verified.

4.3 Cleaning

Clean the support regularly.

- 1) Remove the plastic splint.
- 2) **Recommendation:** Use a laundry bag or net.
- 3) Wash the support at 30 °C in warm water with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

6.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

6.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH