







3R60-PRO, 3R60-PRO=ST, 3R60-PRO=KD, 3R60-PRO=HD





DE	Gebrauchsanweisung	6
EN	Instructions for use	15
FR	Instructions d'utilisation	24
IT	Istruzioni per l'uso	33
ES	Instrucciones de uso	43
PT	Manual de utilização	53
NL	Gebruiksaanwijzing	62
SV	Bruksanvisning	72
DA	Brugsanvisning	80
NO	Bruksanvisning	89
PL	Instrukcja użytkowania	97
HU	Használati utasítás	107
CS	Návod k použití	116
TR	Kullanma talimatı	125
EL	Οδηγίες χρήσης	133
RU	Руководство по применению	144
JA	取扱説明書	154
ZH	使用说明书	162



▲	■
(D) Mindestmenge	Einzelteile
(GB) Minimum order quantity	Single components
(F) Quantité minimum	Pièces à l'unité
(I) Quantità minima	Singoli componenti
(E) Cantidad mínima	Componentes
(P) Quantidade mínima	Componentes
(NL) Minimikvantitet	Onderdelen
(S) Minimum aantal	Separata delar
(DK) Mindste mængde	Komponenter
(N) Minste antall	Separate deler
(PL) Ilość minimalna	Podzespoły
(H) Minimum mennyiség	Alkatrészek
(CZ) Minimální množství pro objednání	Jednotlivé díly
(RO) Cantitate minimă	Componente individuale
(HR) Minimalna količina	Pojedinačni dijelovi
(TR) Minimum miktar	Münferit parçalar
(GR) Ελάχιστη ποσότητα	Μεμονωμένα εξαρτήματα
(RUS) Минимальное количество	Отдельные детали
(J) 複数入パックで発注いただく部品	1個から発注いただける部品
(中国) 最低起订量	单个部件

INFORMATION

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INFORMATION

Last update: 2015-05-07

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1 Single components (Fig. 1)

- (1) 710H10=2X3 Adjustment Wrench
- (2) 4Z59=5.5X6 Extension Stop
- (3) 4G444=3R60-PRO Signature Plate
- (4) 4G432 Upper Protective Knee Cover
- (5) 4G433 Lower Protective Knee Cover
- (6) 4G70 Lamination Anchor with
 - 4 pcs of 501T1=M5×16 Truss Head Screw (a)
 - 4 pcs of 502R1=M5×16 Round Nut with two holes (b)
 - 2 pcs of 506G3=M8X12-V Set Screw (c)
 - 1 pc. of 506G3=M8X10 Set Screw (d)


2 Description

2.1 Indications for use

The **EBS^{PRO}** Knee Joint is to be exclusively used for the prosthetic fitting of amputations of the lower limb. It is a polycentric joint offering a hydraulically controlled, ergonomically balanced stride and hydraulic swing phase control.

2.2 Field of application

Field of application according to the **Otto Bock MOBIS Mobility System**:

Recommended for amputees with **Mobility Grades 2 and 3** (restricted outdoor walkers, non-restricted outdoor walkers).
Approved for a **patient weight of up to 75 kg / 165 lbs.**

The Otto Bock MOBIS is not intended to be used as a guide to obtaining reimbursement for prosthetic components in the USA.

Caution!

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.
- ▶ **Inform the patient.**

2.3 Environmental conditions

Attention!

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water, acids and other liquids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

Please inform your patients.

Allowable environmental conditions

Temperature range for use: -10 °C to +60 °C

Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions

Mechanical vibrations or impacts

Perspiration, urine, fresh water, salt water, acids

Dust, sand, highly hygroscopic particles (e. g. talcum)

2.4 Construction and function (Fig. 2–4)

The amputee will experience considerably greater comfort and stability during stance and swing phase with the proven joint. Two axial forks connect the upper and lower joint section to one another and create a kinematic chain. A pivoting mount (Fig. 2, A) joins the lower joint section to the rear axial fork, which is also connected to the *EBS^{PRO}* shock absorbing unit (Fig. 2, B).

At heel strike, the proximal joint sections swing posteriorly around the lower axes (Fig. 4, D). The *EBS^{PRO}* unit is compressed as the pivoting mount rotates (Fig. 3b and 4E). This pivoting mount serves as a visual control for measuring the efficiency of the *EBS^{PRO}* function. The resistance of the *EBS^{PRO}* unit is continuously adjustable and can be set for the weight and activity level of the patient (see Section 4.3.2.2).

The multiaxial action of the joint is composed of rotary and gliding movements. The center of rotation (instantaneous center of rotation) varies with the flexion position (Fig. 4, F).

The *EBS^{PRO}* unit (EBS = Ergonomically Balanced Stride) is a special feature. With special, hydraulically created resistances the unit at heel strike allows cushioned stance flexion of up to 15° before normal swing flexion is initiated. The unit's progressive cushioning allows optimal control of stance phase flexion in dependence of the walking speed of the prosthesis wearer. When walking slowly, the *EBS^{PRO}* unit allows a higher degree of knee flexion at heel strike. In conjunction with its polycentric kinematics, the joint thus provides optimal stability and comfort. With higher walking speeds stance flexion damping increases, thereby increasingly limiting flexion and enabling the amputee to walk dynamically and in energy efficient manner. **Thanks to this function, walking with the prosthesis becomes much more comfortable and physiological over a large range of walking speeds.**

Another advantage is the very large knee flexion angle, which goes far beyond the usual range. The joint does not feature a flexion stop bumper. A flexion angle of more than 175° is technically pos-

sible, and is limited only by the corresponding socket connectors, socket shapes or by the cosmetic foam cover.

An individually adjustable hydraulic unit located between the axial forks controls swing phase. It prevents excessive heel rise and dampens extension before initial impact. The hydraulic unit's resistances are independently adjustable.

2.5 Different connection systems

The *EBS^{PRO}* Knee Joint is available in 4 different versions. They differ from each other solely by the type of proximal connection used (see front page):

Order number	Version
3R60-PRO	Knee Joint with pyramid connector
3R60-PRO=ST	Version with threaded connector for long residual limbs
3R60-PRO=KD	Version with lamination anchor for knee disarticulation amputees
3R60-PRO=HD	Version with pyramid connector anteriorly angled by 10° for hip disarticulation amputees

Attention!

For the prosthetic fitting of hip disarticulation amputees the specially adapted 3R60-PRO=HD joint version must be used. The proximal joint connector facilitates prosthetic alignment of the anteriorly placed hip joint.

3 Technical data

Order number	3R60-PRO	3R60-PRO=KD	3R60-PRO=ST	3R60-PRO=HD
Proximal connection	Adjustment pyramid (movable)	Lamination anchor	Threaded connector	Adjustment pyramid, angled by 10°
Distal connection	Adjustment pyramid			
Knee flexion angle	175°	145°	125°	175°
Weight	770 g	840 g	750 g	770 g
System height	150 mm	169 mm	165 mm	150 mm
Proximal system height up to alignment reference point	2 mm	21 mm	17 mm	2 mm
Distal system height up to alignment reference point	148 mm			
Max. user weight	75 kg / 165 lbs			
Mobility Grade	2, 3			

4 Handling

4.1 Alignment

The three-dimensional relationship between the socket and the endoskeletal parts affects the dynamic and static functions of the prosthesis. Function of the joint is influenced by the position of the axes. In extension, the instantaneous center of rotation is located above the adjustment pyramid and behind the alignment reference line, thus providing knee stability in the stance phase (Fig. 4). The advantages of the *EBS^{PRO}* Knee Joint can only be made optimal use of in case of a correct alignment.

The optimal residual limb position must be anticipated when positioning the socket connector. Plumb lines in the frontal and sagittal planes drawn from the hip joint's center of rotation and marked during plaster cast taking and trial fitting of the test socket will facilitate correct positioning of the lamination anchor/ socket adapter.

To align the prosthesis please proceed in two steps:

1. Bench alignment using an alignment tool such as 743L200 L.A.S.A.R. Assembly.
2. Static alignment optimization using 743L100 L.A.S.A.R. Posture.

4.1.1 Bench alignment with alignment tool (the following steps refer to Fig 12)

- ① Position the middle of the foot **30 mm anterior** to the alignment reference line.
- ② Set the required heel height of the foot and add 5 mm. Set correct outward rotation of the foot.
- ③ Clamp the knee joint using the appropriate adapter inserts. For bench alignment, the **alignment reference line should run through the upper anterior axis (alignment reference point)**. At that point the pyramid base should be horizontal. Pay attention to the knee ground distance and outward rotation of the knee (adapter inserts provide for a rotation of approx. 5°). Recommended positioning of the alignment reference point: 20 mm above the medial tibial plateau.
- ④ Connect foot and modular knee joint using a tube adapter.
- ⑤ Mark the center of the socket proximal and distal on the lateral side. Draw a line through both marks from socket brim to the distal end of the socket.
- ⑥ Now position the socket such that the alignment reference line passes through the proximal center mark of the socket. Set the socket flexion to somewhere between 3° and 5°; however, the individual situation (e.g. hip joint contractures) must be taken into account and, if necessary, more flexion should be provided. Also pay attention to the ischial tuberosity to ground distance.

Attention!

If the residual limb flexion is not taken into account, the joint will be positioned too far to the front. This will lead to malfunction and premature wear. Use the shiftable adjustment pyramid (4G235=T) of the *EBS^{PRO}* Knee Joint or, in addition, the 4R118 Adapter Plate to optimize the alignment or for later alignment corrections.

Adjusting the movable adapter

The 4G235=T Connection Adapter is designed exclusively for the use with the 3R60-PRO Otto Bock Knee Joint. The adapter allows for displacement of maximally 10 mm in A/P direction, e.g. for subsequent optimization of the prosthetic alignment (Fig. 11a).

To make an adjustment, you must proceed as follows (Fig. 11b):

1. First bring the knee joint into full flexion and remove the cover cap (Important: Keep in mind the position of the cover cap before removing it);
2. Subsequently, loosen the screw just enough that the connection adapter allows for shifting;
3. Then shifting can be made as desired for optimization of the prosthetic alignment;

4. Finally, tighten the screw with 25 Nm and insert the cover cap. It is very important to insert the cover cap in the correct position.

- ⑦ Connect socket and knee using a corresponding adapter (e.g. 4R111, 4R41, 4R55, 4R51 Socket Adapter).

4.1.2 Static alignment optimization using 743L100 L.A.S.A.R. Posture

(the following steps refer to Fig. 13)

The bench alignment can be substantially improved using L.A.S.A.R. Posture. In order to ensure an appropriate stability and at the same time an easy swing phase initiation, proceed as follows:

- ① To make the load line visible, the transfemoral amputee stands on the L.A.S.A.R. Posture with the prosthetic side on the force plate and with the other leg on the height compensation panel. The prosthesis side should be sufficiently loaded (> 35 % of the body weight).
- ② Now adapt the alignment by only **adjusting the plantar flexion** of the foot. The **load line** (laser line) should be approx. **10 mm anterior to the lower front knee axis** (see Fig. 13).
- ③ After step 2, dynamic optimization can take place between parallel bars (see Section 5).

4.2 Combination possibilities

Depending on the patient's functional demands, the following prosthetic feet can be used: Adjust (1M10), 1A30 Greissinger plus, 1D35 Dynamic Motion, 1C30 Trias, 1C40 C-Walk, 1E56 Axion, 1E58 AxionDP or 1E42/43 LuXonMax DP.

4.3 Adjustments and final assembly

4.3.1 Socket fabrication 3R60-PRO=KD

4.3.1.1 Lamination prior to trial fitting

Take a piece of 623T3 Perlon Stockinette twice the length of the plaster model and pull half of it over the isolated socket. Twist the other half of the stockinette on the distal end and reflect it over the socket. In order to absorb strong forces in the area of the modular knee joint, the socket must be gradually reinforced with 616G13 Fiberglass Stockinette. To do this, pull over the first layer to 2/3 of the socket length, tie it off and then reflect it up to half of the socket length. Put 2 layers of 616G12 Carbon-Fiber Cloth on the distal area so as to provide an underlying base with a projecting brim of 3 cm around the 4G70 Adapter, which will be mounted later in the procedure. Pull on two layers of 623T3 Perlon Stockinette. The socket lamination is performed in two stages. In the first stage the socket is laminated up to 2/3 of its length with 617H19 Orthocryl Lamination Resin. In the second stage the proximal part of the socket is laminated with 617H17 Orthocryl Flexible Resin. Once the first casting has hardened, pull on 2 more layers of 623T3 Perlon Stockinette before laminating the proximal part of the socket with 617H17 Orthocryl Flexible Resin.

Attention!

The lamination anchor acts as a flexion stop in the area of the *EBS^{PRO}* pivoting mount (Fig. 10). To provide this stop function, and in order to avoid damage to the knee joint, an even contact surface in the area of the *EBS^{PRO}* pivoting mount must be present after final lamination.

4.3.1.2 Applying the lamination anchor

The correctly fitted lamination anchor must be prepared before it is glued to the socket. To do this, fill the coupling opening with 636K8 Plastaband. Clamp the socket and joint into the alignment apparatus.

Remove soft, flexible or porous materials from beneath the support surface of the lamination anchor. Mix the putty exclusively from 617H21 Orthocryl Sealing Resin and 639A1 Talcum Powder, and glue on the lamination anchor.

Secure with 627B2 Adhesive Tape for trial fitting. Test the stop function (see Section 4.1). If necessary, form the required contact surface with putty. It is also possible to glue a stop bumper of Pedilin onto the outer laminate.

4.3.1.3 Socket completion

Once the trial fitting has been completed, attach the lamination anchor to the socket with truss head screws and joint nuts, and then laminate.

Proceed with further reinforcement as follows:

1. Two layers of 616G12 Carbon Fiber Cloth must be sandwiched between two layers of 623T3 Perlon Stockinette. To do this, pull one layer of 623T3 Perlon Stockinette over the whole socket and tie it on the top circularly, apply the two layers of 616G12 Carbon Fiber Cloth over the arms of the 4G70 Adapter, and then reflect the remaining half of the 623T3 Perlon Stockinette over the carbon fiber layers.
2. Now the socket is to be reinforced once again with 616G13 Fiberglass Stockinette (as described in Section 4.3.1.1).
3. Finally, pull on 2 layers of 623T3 Perlon Stockinette. Proceed with lamination as performed in the first stage of lamination.

Attention!

Deviating from the processing instructions and the recommended materials for the lamination of the anchor can lead to a loosening or breaking of the adapter.

Upon completion of the assembly tighten the adjustment screws with 710D1 Torque Wrench. Torque: 15 Nm/11 lbf.ft/133 lbf.in. Once you have attained the final setting for the prosthesis, use 636K13 Loctite to secure the set screws.

4.3.2 Adjusting the knee joint during trial walking

The user will notice that the *EBS^{PRO}* Knee Joint is different from other prosthetic knees because of the Ergonomically Balanced Stride function. With *EBS^{PRO}*, the user can walk with a natural gait, allowing the knee to flex during stance phase without risk of falling (Fig. 3a). Trial walking should be started with the bench alignment and the factory-set base setting. During trial walking, explain the important, biomechanically advantageous *EBS^{PRO}* function to the patient.

For optimal function, in addition to proper alignment and adjustments, the patient must be taught to use the stance flexion feature. The comfort and stability provided by the knee are most noticeable when walking long distances or on uneven surfaces. Initially, the patient may be unaccustomed to the cushioned knee flexion at heel strike. With proper instruction and experience, the amputee will learn to take considerable advantage of this more natural gait pattern. Therefore, allow the prosthesis wearer to get accustomed to the joint.

Before changing the factory settings, the following guidelines must be observed:

- For changing the degree of stance flexion, first correct the position (see Section 4.3.2.1), then re-adjust the *EBS^{PRO}* unit (see Section 4.3.2.2).

- For swing phase, adjust flexion resistance first and then adjust extension resistance (see Section 4.3.2.3).

4.3.2.1 Stance phase stability by positioning the joint (Fig. 5+6)

Contrary to single-axis knee joints, polycentric knee joints are stable at heel strike. In the *EBS^{PRO}* Knee Joint, stance flexion adds additional stance phase stability. At the factory, the *EBS^{PRO}* unit is adjusted to provide medium resistance. Initiation of swing flexion and thus the stability at heel strike depends on the position of the joint and thus of the location of the instantaneous center of rotation. Tipping the joint in the sagittal plane (changing the angle using the adjustment screws) can be used to determine the location of the center of rotation and thus to vary the degree of stability.

Insufficient stance phase stability (knee joint buckles)	=	instantaneous center of rotation (ICR) too far anterior	➡	Tip the joint back-ward using the adjustment screws	(Fig. 5)
Stance phase stability too great (swing phase can only be initiated with a great deal of effort)	=	instantaneous center of rotation (ICR) too far posterior	➡	Tip the joint forward using the adjustment screws	(Fig. 6)

Any adjustments performed above the joint must be corrected by corresponding counter-adjustments below the joint, in order to maintain the position of the foot.

Attention!

The =KD and =ST connection systems require another way of proceeding during manufacture of the prosthesis. Stance phase stability and initiation of flexion depend on the position of the lamination anchor. Later readjustment of the prosthesis in the sagittal and frontal planes is not possible!

4.3.2.2 Adjusting the stability of the elastic footfall

The degree of elastic stance flexion resistance may be adapted by adjusting the *EBS^{PRO}* unit with the adjustment nut. To do this, insert 710H10=2x3 Adjustment Wrench provided with the knee into the hole (Fig. 7).

Rotation to the left (in - direction)	= decreased resistance	= easier knee flexion	= strong raising of the pivoting mount (Fig. 3b)	= increased stance stability
Rotation to the right (in + direction)	= increased resistance	= less knee flexion	= slight raising of the pivoting mount (Fig. 3a)	= decreased stance stability

Attention!

If the maximally adjustable resistance of the EBS unit is insufficient (knee flexion too easy), possibly tip the joint forward (see Section 4.3.2.1), but be aware that knee stability is reduced by this measure.

4.3.2.3 Adjusting the damping for swing phase control

The hydraulic swing phase control ensures a smooth gait. Movement resistances prevent the prosthetic lower leg from swinging too far when flexing and guarantee dampened extension. Along with

the patient's habits, the foot weight and the length of the lower leg, which have a pendular weight, also have an influence on his or her gait. Upon delivery, the valve screws are set to the lowest resistance for extension damping. For flexion damping, a medium setting has been chosen. The flexion and extension resistances can be independently adjusted using the 710H10=2x3 Adjust-ment Wrench provided.

Adjusting flexion (F) on the left joint side (from the posterior direction of view, Fig. 8).

Rotation of the valve screw of the hydraulic unit to the right (+)	= greater resistance	= less heel rise
Rotation of the valve screw of the hydraulic unit to the left (-)	= less resistance	= more heel rise

Adjusting extension (E) on the right joint side (Fig. 9).

Rotation of the valve screw of the hydraulic unit to the right (+)	= greater resistance	= more extension damping
Rotation of the valve screw of the hydraulic unit to the left (-)	= less resistance	= less extension damping

Attention!

Proceed **with caution and in small increments** when adjusting the resistances, since there is the risk of falling!

When making the final extension damping adjustment, make sure that the leg still completely extends with each step. When adjusting the swing phase control, the influence of a cosmetic foam cover must be compensated for.

4.4 Cosmetic foam cover

Use the 3S107 Foam Cover for the **EBS^{PRO}** Knee Joints. Fabrication of a customized Otto Bock cosmetic foam cover is possible.

Attention!

Do not use talcum powder to eliminate noise in the cosmetic foam cover. Talcum powder reduces the lubrication of the mechanical parts, which may lead to a malfunction and thus increase the risk of failure. Using this medical product after application of talcum powder will render all claims null and void.

Note:

To reduce friction and to eliminate noise, please apply 519L5 Silicone Spray directly on the contact surfaces of the cosmetic foam cover.

4.5 Maintenance instructions

Attention!

Do not use any aggressive cleaning agents since they could cause damage to bearings, seals and plastic parts.

Do not disassemble the knee joint! If you have a problem, please send the knee joint in for service.

**Attention!**

The knee joint must not be lubricated or greased, since this would result in a risk of bearing damage and loss of function.

**Attention! - Please inform your patients!**

The knee joint function may be affected by both environmental and service conditions. To reduce risk for the patient, the knee joint should not be used when noticeable functional changes occur.

Noticeable functional changes may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc. In extreme cases of dynamic permanent loading, the hydraulic damping systems can heat up.

When using the knee joint do not reach into the joint mechanism since there is the risk of injury.

Measure after noticeable functional changes have occurred:

Have the prosthesis examined at a specialist workshop.

Otto Bock recommends readjusting the knee joint's settings once the patient has gotten used to the prosthesis. This period of time varies depending on individual patient characteristics.

Please check the knee joint for wear and proper function at least once a year and make adjustments if necessary. Special attention should be paid to the flexion and extension resistances, bearings, and abnormal noises. Full flexion and extension must be guaranteed.

This component has been tested according to ISO 10328 standard for three million load cycles. Depending on the amputee's activity this corresponds to a duration of use of three to five years.

We recommend carrying out regular yearly safety checks.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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[illegible]

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

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Otto Bock HealthCare Canada
Burlington, Ontario, L7L 5N5
Canada
T +1 289 288-4848
F +1 289 288-4837
infocanada@ottobock.com

Otto Bock HealthCare
Andina Ltda.
Bogotá · Colombia
T +57 1 8619988
F +57 1 8619977
info@ottobock.com.co

Otto Bock de Mexico S.A. de C.V.
C.P. 01180 México, D.F. · Mexico
T +52 55 5575 0290
F +52 55 5575 0234
info@ottobock.com.mx

Otto Bock HealthCare
Minneapolis, MN 55447 · USA
T +1 763 553 9464
F +1 763 519 6153
usa.customerservice
@ottobockus.com

Asia/Pacific

Otto Bock Australia Pty. Ltd.
Baulkham Hills NSW 2153
Australia
T +61 2 8818 2800
F +61 2 8814 4500
healthcare@ottobock.com.au

Beijing Otto Bock
Orthopaedic Industries Co., Ltd.
Beijing, 100015, P.R. China
T +8610 8598 6880
F +8610 8598 0040
news-service@ottobock.com.cn

Otto Bock Asia Pacific Ltd.
Kowloon, Hong Kong · China
T +852 2598 9772
F +852 2598 7886
info@ottobock.com.hk

Otto Bock HealthCare
India Pvt. Ltd.
Mumbai 400 021 · India
T +91 22 2274 5500 / 5501 / 5502
information@indiaottobock.com

Otto Bock Japan K. K.
Tokyo, 108-0023 · Japan
T +81 3 3798-2111
F +81 3 3798-2112
ottobock@ottobock.co.jp

Otto Bock Korea HealthCare Inc.
137-070 Seoul · Korea
T +82 2 597-3831
F +82 2 577-3828
info@ottobockkorea.com

Otto Bock
South East Asia Co., Ltd.
Bangkok 10900 · Thailand
T +66 2 930 3030
F +66 2 930 3311
obsea@ottobock.co.th

Other countries

Otto Bock HealthCare GmbH
37115 Duderstadt · Germany
T +49 5527 848-1590
F +49 5527 848-1676
reha-export@ottobock.de

ISO 10328–P5-75 kg*)

*) Body mass limit not be exceeded !
For further details see manufacturer's
written instructions on intended use !

Otto Bock HealthCare GmbH
Max-Näder-Straße 15 · 37115 Duderstadt / Germany
T +49 5527 848-0 · F +49 5527 72330
healthcare@ottobock.de · www.ottobock.com

Ottobock has a certified Quality Management System in accordance with ISO 13485.