
50K30 Xeleton

DE Gebrauchsanweisung	6
EN Instructions for Use	11
FR Instructions d'utilisation	16
IT Istruzioni per l'uso	21
ES Instrucciones de uso	26
PT Manual de utilização	31
NL Gebruiksaanwijzing	36
SV Bruksanvisning	42
DA Brugsanvisning	46
NO Bruksanvisning	51
PL Instrukcja użytkowania	56
HU Használati utasítás	61
CS Návod k použití	66
SK Návod na používanie	71
RU Руководство по применению	76


18

19

	Seite/ Side	Größe/ Size	Körpergröße/ Body height		Umfang/Circumference *	
			<1,83 m (6 ft)	>1,83 m (6 ft)	cm	inch
	L/R	XS	2	1	33 – 39	13 – 15.4
	L/R	S	2	1	39 – 47	15.4 – 18.5
	L/R	M	2	1	47 – 53	18.5 – 20.9
	L/R	L	2	1	53 – 60	20.9 – 23.6
	L/R	XL	2	1	60 – 67	23.6 – 26.4
	L/R	XXL	2	1	67 – 75	26.4 – 29.5
	L/R	XXXL	2	1	75 – 81	29.5 – 31.9

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INFORMATION

Last update: 2013-01-17

- Please read this document carefully.
- Follow the safety instructions.

These instructions for use provide important information on the fitting and application of the 50K30 Xeleton Knee Joint Orthosis.

2 Product Description

2.1 Product Overview (see Figure 1)

Item	Component	Item	Component
A	Joint cover	D	Pull loop
B	Pad retainer	E	Tibia pad
C	Technogel® condyle pad	F	Screw connection

3 Application

3.1 Application

The orthosis is intended **exclusively** for orthotic fittings of the knee and **exclusively** for contact with intact skin.

The orthosis must be used according to the indications and following the instructions of a physician.

3.2 Indications

- After ruptures of the anterior cruciate ligament (ACL) and posterior cruciate ligament (PCL), collateral ligament injuries (MCL, LCL) or after capsular ligament surgery
- Conservative pre and post-operative joint stabilization (e.g. after meniscus suturing or meniscus replacement) with limited range of motion
- Conservative therapy of chronic complex knee joint instabilities

Indications must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

Not known.

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries; inflammation; prominent scars that are swollen; reddening and hyperthermia of the fitted limb; pronounced varicose veins, especially with return flow impairment; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support/orthosis will be applied; sensory and circulatory disorders in the area of the legs, e. g. in case of diabetic neuropathy.

3.4 Effects

The orthosis serves to stabilise the knee joint.

4 Safety

4.1 Explanation of Warning Symbols

 **CAUTION** Warning regarding possible risks of accident or injury

 **NOTICE** Warning regarding possible technical damage

4.2 General Safety Instructions

 **CAUTION**

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety instructions

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

 **CAUTION**

Contact with heat, embers, or fire

Injuries (e.g., burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers, and other sources of heat.
- ▶ **Inform the patient.**

 **CAUTION**

Breakage of the orthosis frame or joint due to improper deformation

Injuries caused by sharp edges, loss of function

- ▶ Avoid adapting the orthosis repeatedly.
- ▶ Do not make any improper changes to the orthosis.
- ▶ **Inform your patient.**

 **CAUTION**

Reuse on other patients and improper cleaning

Skin irritation, formation of eczema, or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

NOTICE

Contact with salves, lotions, or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to salves, lotions, or other products that contain oils or acids.
- ▶ Inform the patient.

5 Handling

INFORMATION

- ▶ The duration of daily wear and overall period of use are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the physician treating you.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any unusual changes occur (e.g., worsening of the complaint).

5.1 Selecting the Size

1. Measure the circumference of the patient's thigh **15 cm (6")** above the center of the patella and determine the orthosis size (see size table).
2. Determine the length of the orthosis based on the patient's body size (see size table). Ottobock recommends the long version of the orthosis for a body size over **1.83 m (6 ft)**.

5.2 Fitting and Application

5.2.1 Adjusting the Range of Motion (Flexion / Extension)

INFORMATION

The joint stops (medial and lateral) must have the same degree settings.

Available extension stops (square grip):	6°, 10°, 20°, 30°, 45°
Available flexion stops (round grip):	6°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

- 1) Remove the condyle pads from the pad retainers (see Fig. 5).
- 2) Remove the adjustment screws from the pad retainers (see Fig. 6).
- 3) Remove the pad retainers (see Fig. 1, item B) and the joint covers (see Fig. 1, item A); see Fig. 7.
- 4) Insert the desired extension and flexion stops (see Fig. 8).
→ The stops will engage with slight resistance.
- 5) Insert the joint covers with the profile edge facing forwards.
- 6) Attach the pad retainers (see Fig. 7).
- 7) Insert the adjustment screws in the pad retainers and tighten them (see Fig. 6).
- 8) Attach the condyle pads to the pad retainers with hook-and-loop material (see Fig. 5).
- 9) Move the orthosis joints, and check the position of the stops.

5.2.2 Adaptation to the Knee Width

Spacer pads (thickness = 5 mm and 8 mm) are included in the scope of delivery. The spacer pads can be mounted between the pad retainers (see Fig. 1, item B) and the condyle pads (see Fig. 1, item C) using hook-and-loop fasteners (see Fig. 5). This adjustment can be carried out independently on the medial and lateral joints.

5.2.3 Applying the ACL Version

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.
- ▶ **Inform the patient.**

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear, and damage prior to each use.
- ▶ Advise the patient that the product may not be used if the product or part of the product exhibits signs of wear (e.g., cracks, deformation, poor fit) or damage.

- > All straps are opened.
- > The patient is sitting on the front edge of a chair. The knee is angled at 45°.
- 1) Position the orthosis with the middles of the joints approximately even with the top of the patella (see Fig. 9).
 - Although the orthosis starts out being a little too high, it will quickly settle into the correct position.
- 2) Fasten the lower front strap for the lower leg (see Fig. 10).
- 3) Fasten the upper rear thigh strap (see Fig. 11).
- 4) Keep the knee bent, and push the two joints towards the back of the leg (see Fig. 12).
- 5) Fasten the lower thigh strap (see Fig. 13). Pull the joints forward slightly, and position them so they are centered relative to the knee joint. Ensure that the joints are not pulled ahead of the lateral center line of the leg (see Fig. 14).
- 6) Fasten the upper front strap for the lower leg (see Fig. 15).
- 7) Fasten the rear strap for the lower leg (see Fig. 16).
- 8) Tighten the orthosis straps so that the orthosis fits firmly but is comfortable to wear.
- 9) After **15 min** of active movement (e.g., slight knee bends): slightly retighten the straps in the specified sequence (see Fig. 2, straps 1-5).

5.2.4 Conversion to CI Version

- > The orthosis is in the as-delivered condition (ACL version).
- > The CI accessory set is on hand.
- 1) Loosen and remove the screw connection on the lower thigh strap (strap 3); see Fig. 17.
- 2) Attach the pull loop on the lower thigh strap and the additional pull loop (see Fig. 1, item D; CI accessory kit) with the longer screw (see Fig. 18).
- 3) Mount the additional strap (strap 6, CI accessory kit) as the front thigh strap.
- 4) Position the orthosis correctly (see the section "Applying the ACL Version").
- 5) Fasten the straps in the specified sequence (see Fig. 3, straps 1-6).

5.2.5 Conversion to PCL Version

- > The orthosis is in the as-delivered condition (ACL version).
- > The PCL accessory kit is on hand.
- 1) Remove the lower thigh strap (strap 3), the rear strap for the lower leg (strap 5), and the upper front strap for the lower leg (strap 4).
- 2) Turn the pull loop on the lower thigh strap (see Fig. 1, item D) to the front.
- 3) Use strap 4 without tibia pad with the pad on strap 3 as the front thigh strap.

- 4) Loosen the screw connection of the upper front strap for the lower leg, and remove the front pull loop (see Fig. 19).
- 5) Mount the pull loop on the rear strap for the lower leg with the short screw (PCL accessory kit).
- 6) Use strap 3 with the pad on strap 5 as the rear strap for the lower leg.
- 7) Position the orthosis correctly (see the section "Applying the ACL Version").
- 8) Fasten the straps in the specified sequence (see Fig. 4, straps 1-4).

5.2.6 Replacing the Straps and Pads

Straps:

- The straps are equipped with opposing hook-and-loop closures and are easy to replace.

Strap and frame pads:

- The strap and frame pads are attached with hook-and-loop closures and are easy to replace.

5.2.7 Use of undergarments and covers

Optional undergarments and covers are available for the orthosis. An undergarment is used as a layer between the orthosis and the skin, therefore protecting against skin damage. A cover is pulled over the orthosis, protecting it against damage. This also reduces the risk of injury in sports with the possibility of massive body contact (e.g. collisions).

The following products can be used with the orthosis (the size is specified in the reference number):

Designation	Reference number
UNDERGARMENT FOR ACL KNEE ORTHOSIS	29K193=M/L/XL
UNDERGARMENT FOR ACL KNEE ORTHOSIS	29K193=XS/S
UNDERGARMENT FOR ACL KNEE ORTHOSIS	29K193=XXL/XXXL
TEXTILE COVER FOR ACL KNEE ORTHOSIS	29K194=M/L/XL
TEXTILE COVER FOR ACL KNEE ORTHOSIS	29K194=XS/S
TEXTILE COVER FOR ACL KNEE ORTHOSIS	29K194=XXL/XXXL

5.3 Delivery

Check for proper fit of the orthosis on the leg at the time of delivery. Ensure that the centers of the patella and the joints are approximately lined up with each other.

5.4 Cleaning

INFORMATION

- Ensure proper alignment when inserting the tibia pad. The wide side of the pad has to face up.

- 1) Remove the straps and pads.
- 2) Hand-wash the straps and pads in warm water at 30°C (86°F) using neutral detergent.
- 3) Allow to air-dry. Do not expose to direct heat (e.g., from the sun, an oven or radiator).
- 4) If the orthosis comes into contact with salt water or dirt, rinse with clear water and allow to air-dry.

6 Legal Information

6.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

6.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This product has been classified as a class I device according to the classification criteria outlined in appendix IX of the directive. The manufacturer therefore drew up the declaration of conformity under its own responsibility according to annex VII of the directive.

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Otto Bock HealthCare GmbH