


8065

DE	Genu Immobil Vario Knieorthese	2
EN	Genu Immobil Vario Knee Orthosis	4
FR	Orthèse de genou Genu Immobil Vario	6
IT	Ortesi di ginocchio Genu Immobil Vario	8
ES	Rodillera Genu Immobil Vario	10
NL	Genu Immobil Vario knieorthese	12
DA	Genu Immobil Vario knæortose	14
PL	Orteza stawu kolanowego Genu Immobil Vario	16
HU	Genu Immobil Vario térdortézis	18
CS	Kolenní ortéza Genu Immobil Vario	20
SK	Ortéza kolena Genu Immobil Vario	22
EL	Νάρθηκας γόνατος Genu Immobil Vario	24
RU	Коленный ортез Genu Immobil Vario	26
ZH	Genu固定可调膝关节矫形器	28

Material	Al; Baumwolle/Cotton; PA; PES; PP; PU; Schmelzkleber/Hotmelt adhesive; Stahl/Steel (Nickel Legierung/nickel alloy)
-----------------	--

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Deutsch

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1 Application

English

The 8065 Genu Immobil Vario is to be used **solely** for the orthotic fitting of the knee joint.

2 Indications and Effects

Post-traumatic and post-operative immobilization and mobilization of the knee joint through limitation of movement. Indication must be determined by the physician.

2.1 Contraindications

2.1.1 Absolute Contraindications

None.

2.1.2 Relative Contraindications

In case of the following indications, consultation with a physician is required: Skin diseases/injuries, inflammatory symptoms, prominent scars on body parts to which the support is to be applied that are swollen, reddened and feel excessively warm; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support is to be applied; sensory and circulatory disorders of the leg.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warnings regarding possible risks of accident or injury.

 **NOTICE** Warnings regarding possible technical damage.

3.2 General Safety Instructions

CAUTION

The orthosis may be over-tightened

Risk of local pressure and constriction of blood vessels and nerves.

► Verify the correct positioning and fit of the orthosis.

CAUTION

Use on more than one patient

Skin irritation, formation of eczema or infections

► Only use the orthosis on one patient.

NOTICE

Incorrect environmental conditions

The orthosis is not flameproof.

- ▶ Do not expose the orthosis to open flame or other heat sources.

NOTICE

The lifespan of the materials may be reduced

- ▶ Avoid contact with substances containing oils and acids, salves and lotions.

4 Trial Fitting and Delivery of the Orthosis

Upon delivery of the orthosis to the patient, check for proper fit on the leg. The orthotic joints shall be positioned at mid-patella height. Otherwise, adjust the joint on the limb as appropriate (in the compromise center of rotation in accordance with the indications by Niethard).

5 Adjusting the Joints

In order to set the joint free, you have to push upwards the releasing button; you must use the plugs to obtain the wished movement radius. Push therefore down the releasing button to fix the chosen position (fig. 1).

⚠ CAUTION

Improper adjustment of the orthotic joints

Risk of injury and damage to the orthosis

- ▶ Make sure that the adjustment of both orthotic joints is equal.

6 Fitting and Donning the Orthosis

Lay the femoral and tibial pads around the leg and frontally close the hook & loop closures (fig. 2).

Position the lateral joint at mid-patella height and attach it with the hook strip (on the joint bar's back) to the pads (fig. 3). Use the same procedure to attach the medial joint.

⚠ CAUTION

Improper changes of joint settings

- ▶ Adjustments of the joint bars may only be made by qualified personnel.
- ▶ Adjust the orthosis so that the joint bars are parallel to each other and that the joints are at the same height.

Pull the straps through the loops, tighten the straps and attach them with the hook & loop closures (fig. 4). For doffing the orthosis, just open the straps and pads without detaching the joint bars. This will avoid repeated adjustment.

7 Care Instructions

Hand wash the pads in lukewarm water with a standard mild detergent. Thoroughly rinse and let air dry. Do not tumble dry. No chemical cleaning.

8 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

9 CE Conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of

the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH