Otto Bock®

QUALITY FOR LIFE

4165 Manu Neurexa

(D)	Handgelenkorthese	3
(GB)	Wrist Orthosis	7
(F)	Orthèse de poignet	10
(I)	Ortesi di polso	14
(E)	Órtesis de muñeca	18
(P)	Ortótise para punho	22
(NL)	Polsorthese	26
(S)	Handortos	30
(DK)	Håndledsortose	34
(N)	Håndleddortose	37
(FIN)	Ranneortoosi	40
(PL)	Orteza stawu nadgarstka	44
(H)	Csuklóórtézis	48
(CZ)	Ortéma zápěstí	52
(GR)	Νάρθηκας καρπού	56
(J)	手首装具	60

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⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

INFORMATION

Please read these instructions for use carefully and pay special attention to the safety information!

1 Description

1.1 *Intended Use*

Manu Neurexa is to be used **exclusively** for orthotic treatment of the wrist and is **exclusively** suitable for contact with healthy skin.

1.2 *Indications and Effects*

1.2.1 Indications

Paresis of the forearm muscles resulting in weakness of the wrist e.g. after stroke, cervical disk herniation, cervical plexus injury, peripheral nerve damage. Indication must be determined by the physician.

1.2.2 Effects

Manu Neurexa stabilizes the wrist and hand. It supports the sensorimotor control of the hand and facilitates occupational training, hand function and physical activity.

1.3 *Contraindication*

1.3.1 Absolute contraindication

Severe spasticity in the arm and wrist.

1.3.2 Relative contraindication

In case of the following indications, consultation with a physician is required: Skin diseases / -injuries, inflammatory symptoms, prominent scars that are swollen, reddened and feel excessively warm when touched; Lymphatic flow disorders, including unclear soft tissue swelling distant to the body area

to which the appliance will be applied; Sensory and circulatory disorders of the arm and hand.

1.4 Safety Instructions

INFORMATION

The patient is to be instructed in the proper use / care of the product.

The initial fitting and application of the product must be carried out by trained, qualified personnel only. The daily duration of use and period of application are dependent on medical indication by the physician.

⚠ CAUTION

Risk of injury as a result of improper use. The product is designed for use on one patient only. Parts to be fitted and those parts that come directly into contact with the skin can cause functional and hygienic risks if the orthosis is used by another person.

An orthosis / support applied too tightly to the body can cause local pressure and, in some cases, even restrict adjacent vessels or nerves. Do not apply the product too tightly. Consult a physician immediately if you experience unusual changes (such as increase in pain).

Improper changes to the product are not permitted.

NOTICE

Risk of damage due to incorrect environmental conditions. This product is not flame-resistant. Keep the product away from flames or other heat sources.

The product should not come into contact with grease or acidic agents, unguents and lotions, as this can reduce the useful life – especially of pads.

2 Handling

2.1 Size Selection of the product

The size selection of the orthosis is based on the wrist circumference.

2.2 Adapting and Application of the product

1. Open the orthosis.

2. Put the wrist orthosis on a table, plastic splint side down (Fig. 1).
3. Place the wrist into the orthosis (Fig. 1).
4. Close the velcro, start with the middle tape, after that the two others (Fig. 2). If necessary correct the middle tape.
5. Make sure the orthosis is firmly fitted without too much pressure.
6. If necessary, adaptation of the pre shaped plastic splint can be done by the fitter / orthotist. This is performed by heating it in boiling water and then adjustment is possible.

2.3 Instruction for use and care

Material: PCM®, plastic splint (Polyamide), steel splints, Polyester, velcro.

Cleaning: Machine wash at 40 °C / 104 °F with a standard mild detergent. Wash out thoroughly and let air-dry. Take out the plastic splint and close all velcros before machine wash. We recommend using a washing bag. Note: Residues of cleaning agents may cause skin irritations and wear of material. To extend the products life-time: close all velcro when not using the product.

Disposal after usage: The material used in the product is recyclable and should be disposed of accordingly.

3 Further Usage Restrictions

The product is designed for use on one patient only. Re-use of the product is not permitted. The daily duration of use and period of application are dependent on the medical indication.

4 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

5 CE Conformity

This product meets the requirements of the 93/42/EEC guidelines for medical devices. This product has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH