
8165 Genu Neurexa

DE Gebrauchsanweisung	3
EN Instructions for use	6
FR Instructions d'utilisation	9
IT Istruzioni per l'uso	12
ES Instrucciones de uso	15
PT Manual de utilização	18
NL Gebruiksaanwijzing	21
SV Bruksanvisning	24
DA Brugsanvisning	27
NO Bruksanvisning	30
FI Käyttöohje	33
PL Instrukcja użytkowania	36
HU Használati utasítás	39
CS Návod k použití	42
EL Οδηγίες χρήσης	45
JA 取扱説明書	48

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Explanation of Symbols

English

 CAUTION	Warnings regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.
 INFORMATION	Information regarding operation. Information for service personnel.

INFORMATION

Last update: 2015-02-02

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1 Description

1.1 Intended use

Genu Neurexa is to be used **exclusively** for orthotic treatment of the knee and is exclusively suitable for contact with healthy skin.

1.2 Indications and effects

1.2.1 Indications

Paresis of leg muscles resulting in knee hyperextension e.g. after stroke, lumbar disc herniation, peripheral nerve damage, etc.

Indication must be determined by the physician.

1.2.2 Effects

Genu Neurexa stabilizes the knee in the sagittal and coronal plane and prevents hyperextension and collateral (mediolateral) instability.

It supports the sensorimotor control of the knee and facilitates gait training and physical activity.

1.3 Contraindication

1.3.1 Absolute contraindications

None.

1.3.2 Relative contraindications

In case of the following indications, consultation with a physician is required:

Skin diseases/ -injuries, inflammatory symptoms, prominent scars that are swollen, reddened and feel excessively warm when touched; Lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the appliance will be applied; Sensory and circulatory disorders of the leg.

1.4 Safety instructions

INFORMATION

The patient is to be instructed in the proper use / care of the product.

The initial fitting and application of the product must be carried out by trained, qualified personnel only. The daily duration of use and period of application are dependent on medical indication by the physician.

⚠ CAUTION

Risk of injury as a result of improper use. The product is designed for use on one patient only. Parts to be fitted and those parts that come directly into contact with the skin can cause functional and hygienic risks if the orthosis is used by another person.

An orthosis / support applied too tightly to the body can cause local pressure and, in some cases, even restrict adjacent vessels or nerves. Do not apply the product too tightly. Consult a physician immediately if you experience unusual changes (such as increase in pain).

Improper changes to the product are not permitted.

NOTICE

Risk of damage due to incorrect environmental conditions. The product is not flame-resistant. It must be kept away from flames or other heat sources.

The product should not come into contact with grease or acidic agents, unguents and lotions, as this can reduce the useful life – especially of pads.

2 Handling

2.1 Size selection of the product

The size selection of the orthosis is based on the knee circumference most important 15 cm below midpatella, second hand 15 cm above midpatella.

2.2 Adapting and application of the product

2.2.1 Adapting

The splints are individually adjusted by trained personnel.

1. Adjust the lateral splints.

An orthotist, orthopaedic technician or therapist should adapt the aluminium splints for an individual optimal fit. The splints should be adapted to the curvature of the leg.

INFORMATION

The splint should not be bent less than **3 mm** from the plastic part in the midsection of the splint. This could otherwise reduce the stability of the splint.

2. Adjust desired extension stop and / or flexion stop of the aluminium splints.

Graduated marking on the plastic stop indicates the angle at which the orthosis stops the leg in flexion or extension (functional angle; Fig. 2).

2.2.2 Application

1. Apply the knee orthosis so that the kneecap comes into the centre of the hole at the front. First close the part of the protection which is under the kneecap and then the upper part. Use the materials elasticity and pull up the protection's edges from both directions (Fig. 3).
2. Pull the four straps in the front and close them firmly (Fig. 4).
3. To adjust the extension stop, loosen the crossed velcro fastening in the knee fold. Adjust to the desired length and fasten the velcro (Fig. 5). For maximum extension stop, the knee joint should be bent when you tighten the cross.

2.3 Instruction for use and care

Material: Outlast®, aluminium splints, polyester, velcro.

Cleaning: Machine wash at 40 °C / 104 °F with a standard mild detergent. Wash out thoroughly and let air-dry. Take out the splints and close all velcros before machine wash. We recommend using a washing bag.

Note: Residues of cleaning agents may cause skin irritations and wear of material.

To extend the products life-time: close all velcro when not using the product.

Disposal after usage: The material used in the product is recyclable and should be disposed of accordingly.

3 Further usage restriction

The product is designed for use on one patient only. Re-use of the product is not permitted.

The daily duration of use and period of application are dependent on the medical indication.

4 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

4.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

4.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

4.3 Local legal information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

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Otto Bock HealthCare GmbH