Otto Bock®

QUALITY FOR LIFE

50P12

(D)	Manu Forsa Volar Handgelenkorthese	4
(GB)	Manu Forsa Volar Wrist Orthosis	6
(F)	Orthèse de poignet Manu Forsa Volar	9
(I)	Ortesi di polso Manu Forsa Volar	12
(E)	Órtesis de muñeca Manu Forsa Volar	14
(P)	Ortótise para a articulação da mão Manu Forsa Volar	17
(NL)	Manu Forsa Volar polsorthese	20
(S)	Manu Forsa Volar Handledsortos	22
(DK)	Manu Forsa Volar håndledsortose	25
(N)	Manu Forsa Volar håndleddortose	27
(FIN)	Manu Forsa Volar ranneortoosi	30
(PL)	Orteza nadgarstka Manu Forsa Volar	32
(H)	Manu Forsa Volar csuklóortézis	35
(CZ)	Zápěstní ortéza Manu Forsa Volar	38
(HR)	Ortoza ručnog zgloba Manu Forsa Volar	40
(TR)	Manu Forsa Volar Bilek ortezi	43
(GR)	Νάρθηκας καρπού Manu Forsa Volar	45
(RUS)	Ортез для лучезапястного сустава Manu Forsa Volar	48
(J)	Manu Forsa Volar 手関節装具	51
(中国)	Manu Forsa Volar 腕掌矫形器	53

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1 Intended Use

English

The 50P12 Manu Forsa Volar Wrist Orthosis is to be used exclusively for orthotic treatment of the wrist in accordance with the indications. It is only suited for contact with healthy skin.

2 Indication and Effects

Acute and chronic irritations in the wrist joint (e.g. carpal tunnel syndrome), posttraumatic/post-operative conditions that require the immobilization of the wrist joint, or partial paralyses.

Indication must be determined by the physician. The 50P12 Manu Forsa Volar Wrist Orthosis ensures that the wrist will be largely immobilized in a functional position.

3 Size Selection of the Orthosis

The 50P12 Manu Forsa Volar Wrist Orthosis is delivered in two sizes (for the left and the right side).

The size selection of the orthosis is based on the hand width:

	Size	Hand Width (A)
	S	up to 8.5 cm / 3.35 inch
	M	above 8.5 cm / 3.35 inch

4 Material and Construction

The wrist orthosis is made of skin friendly velour material (polyamide) with a flexible plastic (PE) reinforced insert and a stable metal splint (Nirosta), a velcroed pad covered with velour (polyamide / Lycra®)¹ as well as hook and loop straps (polyamide, fig. 1).

5 Adapting and Applying the Orthosis

The initial fitting and application of the Manu Forsa Volar Wrist Orthosis 50P12 may only be carried out only by qualified personnel.

The Manu Forsa Volar wrist orthosis is delivered in a closed position and ready for a trial fitting.

The metal splint has been anatomically preshaped and is individualized by the orthopedic technician during the adaptation process. If required, the metal splint can be pulled out of the orthosis by opening the appropriate hook and loop closure (Fig. 2+3). If this is done, the appropriate hook and loop closure must be used to close the pocket again after reinserting the metal splint.

The hook and loop are used to secure the pad to the orthosis, which facilitates individual adjustment. The pad is used for the anatomical positioning of the hand and also has a padding effect.

In the forearm area, the posterior hook and loop strap can be adjusted by using the integrated slots and thus individually adapted to the circumference of the forearm.

The middle strap is attached to the orthosis by using the hook and loop strap and is positioned between the ulnar head and the carpal bone (Fig. 7).

The anterior strap is also attached to the orthosis with the hook and loop strap and is set during the adjustment process.

To apply the orthosis, the middle strap is opened and the posterior strap loosened. Then slide the patient's hand into the orthosis from the back (Fig. 4). Afterwards, tighten and fasten first the posterior strap and then the middle strap (Fig. 5+6).

If necessary, re-adjust the anterior strap.

The firm application of all straps, without constriction, improves the fixation effect.

When handing over the orthosis to the patient, the orthosis must be tested for the correct fit on the hand. The patient is also to be instructed in the operation and use of the 50P12 Manu Forsa Volar Wrist Orthosis when it is handed over.

6 Side Effects

Attention!

All aids - supports and orthoses - applied too tightly to the surface of the body may result in local pressure spots or, in some cases, even restrict the adjacent vessels or nerves.

¹ Lycra is a trademark of DuPont

7 Contraindications

Attention!

In the case of the following clinical pictures, the aid should not be applied or worn without first consulting a physician:

1. Skin diseases and injuries to body parts to be fitted, especially in cases of inflammatory symptoms. This also applies to prominent scars that are swollen, reddened and feel excessively warm.
2. Lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the aid will be applied.
3. Sensory disorders - in the arm/ hand

8 Instructions for Use and Care

The orthosis must be used in strict conformity with the indication. The patient's acceptance and cooperation play an important role. It is essential to explain to the patient the proper use of the orthosis! The physician's instructions must be strictly observed.

The daily duration of use and total period of application must be determined by the attending physician.

Patients should be instructed to consult a physician immediately if they experience unusual changes (such as an increase in pain)

Attention!

The orthosis is not flame-resistant. Patients should be instructed to keep the product away from flames or other heat sources.

The orthosis is to be checked carefully for proper function and for possible wear or damage every time it is used. If a hook and loop closure or any other part of the orthosis show any signs of wear, the orthosis should not be used any longer.

The product should not come into contact with grease, acidic agents, ointments or lotions.

Improper changes to the product are not permitted.

Wash by hand in warm water (30° C / 85 °F) with a standard mild detergent. We recommend closing the hook & loop closures before washing in order to prevent any damage to them. To avoid causing skin irritations, rinse the orthosis thoroughly that all soap residues are removed. Do not expose to the sun and do not use oven or radiator heat for drying. Let the orthosis air dry.

9 Usage Restriction Clause

The 50P12 Manu Forsa Volar Wrist Orthosis has been designed for use on only **one** patient. The daily duration of use and period of application are dependent on the medical indication. This orthosis is only suited for contact with healthy skin

Attention!

Parts to be fitted and those parts that come directly into contact with the skin can cause functional and hygienic risks if the orthosis is used by another person.

10 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

11 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH