




50K13N Genu Arexa


DE Gebrauchsanweisung	5
EN Instructions for Use	10
FR Instructions d'utilisation	15
IT Istruzioni per l'uso	20
ES Instrucciones de uso	25
PT Manual de utilização	31
NL Gebruiksaanwijzing	36
SV Bruksanvisning	41
DA Brugsanvisning	46
PL Instrukcja użytkowania	51
HU Használati utasítás	57
CS Návod k použití	62
SK Návod na používanie	67
BG Инструкция за употреба	72
TR Kullanma talimatı	77
EL Οδηγίες χρήσης	82
RU Руководство по применению	87
JA 取扱説明書	93
ZH 使用说明书	98

	Seite/ Side	Größe/ Size	Umfang/Circumference			
			A (cm)	A (inch)	B (cm)	B (inch)
	L/R	S	32–37	12,6–14,6	39–46	15,4–18,1
	L/R	M	37–41	14,6–16,1	46–53	18,1–20,9
	L/R	L	41–46	16,1–18,1	53–61	20,9–24,0
	L/R	XL	46–50	18,1–19,7	61–69	24,0–27,2
	L/R	XXL	50–56	19,7–22,1	69–79	27,2–31,1

Material	AL, SS, klettbarer Velours aus PA/velastic velour of PA, PU-Schaum/PU-foam, PE-Schaum/PE-foam, PA, POM, TPU, PP
-----------------	---

1 Vorwort

INFORMATION

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INFORMATION

Date of the last update: 2013-02-13

- Please read this document carefully.
- Follow the safety instructions.

These instructions for use provide important information on the fitting and application of the 50K13N Genu Arexa knee joint orthosis.

2 Product Description

2.1 Product Overview (see Fig. 1)

No.	Component	No.	Component
1	Strap Guide Loop	5	Technogel® Condyle Pad
2	Pad Button	6	Pad Retainer
3	Clip Rivet	7	Extension Stop (square grip)
4	Tibia Pad	8	Flexion Stop (round grip)

3 Application

3.1 Intended Use

The orthosis is intended **exclusively** for orthotic fittings of the knee and **exclusively** for contact with intact skin.

The orthosis must be used according to the indications and following the instructions of a physician.

3.2 Indication

- After ruptures of the anterior cruciate ligament (ACL) and posterior cruciate ligament (PCL), collateral ligament injuries (MCL, LCL) or after capsular ligament surgery
- Conservative pre and post-operative joint stabilisation (e.g. after meniscus suturing or meniscus replacement) with limited range of motion
- Conservative therapy of chronic complex knee joint instabilities

Indications must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

Not known.

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; pronounced varicose veins, especially with return flow impairment, lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support will be applied; sensory and circulatory disorders in the area of the legs, e.g. in case of diabetic neuropathy.

3.4 Effects

The orthosis serves to stabilise the knee joint.

4 Safety

4.1 Explanation of Warning Symbols

 **CAUTION** Warning regarding possible risks of accident or injury

 **NOTICE** Warning regarding possible technical damage

4.2 General Safety Instructions

 **CAUTION**

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

 **CAUTION**

Contact with heat, embers or fire

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.
- ▶ **Inform the patient.**

 **CAUTION**

Reuse on other patients and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

 **INFORMATION**

- ▶ Observe the information in the section "Information on Re-Use".

CAUTION

Improper moulding or application

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.
- ▶ **Inform the patient.**

NOTICE

Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.
- ▶ **Inform the patient.**

5 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the physician treating you.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

5.1 Size Selection

- 1) Measure the circumference of the knee joint (circumference A).
- 2) Measure the circumference of the thigh (circumference B) 15 cm above the centre of the patella.
- 3) Determine the size of the orthosis (see size table).

5.2 Fitting and Application

5.2.1 Adjusting the Range of Motion (Flexion/Extension)

INFORMATION

The joint stops (medial and lateral) must have the same degree settings.

Available extension stops (square grip):	6°, 10°, 20°, 30°, 45°
Available flexion stops (round grip):	6°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

- 1) Remove the pad retainer (and the spacer rings if applicable) (see Fig. 4).
- 2) Insert the desired extension and flexion stops into the joint (see Fig. 5).
→ The stops will engage with slight resistance.
- 3) **Optional:** Snap the required number of spacer rings onto the joint (see the section "Adaptation to the Knee Width").
- 4) Attach the pad retainer.

5.2.2 Adaptation to the Knee Width

The scope of delivery includes spacer rings (thickness = 5 mm).

This adjustment can be carried out independently on the medial and lateral joints.

- 1) Loosen the pad retainer from the joint.
- 2) Snap the required number of spacer rings onto the joint.
- 3) Attach the pad retainer on the spacer rings.

5.2.3 Adjusting the Strap Height

The height of the front/rear top strap for the lower leg can be adjusted for knee joint relief in case of post-operative therapy and/or to adapt the position of the straps to the anatomical shape of the leg.

- 1) Remove the clip rivets from the front/rear, top strap for the lower leg.
- 2) Move the strap.
- 3) Attach the strap with new clip rivets.

5.2.4 Applying the ACL Version

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- Ensure that the orthosis is applied properly and fits correctly.
- **Inform the patient.**

NOTICE

Use of a worn or damaged product

Limited effectiveness

- Instruct the patient to inspect the product for proper functionality, wear and damage prior to each use.
- Advise the patient that the product must not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

- > The orthosis is in the delivery condition (ACL/CI version).
- > All straps are opened.
- > The patient is sitting on the front edge of a chair. The knee is angled at 45°.
- 1) Position the orthosis with the middle of the joints approximately even with the top of the patella (see Fig. 7).
 - Although this puts the orthosis higher than its final position, it settles into the correct position on its own.
- 2) Fasten the rear, top strap for the lower leg (see Fig. 8).
- 3) Fasten the rear, bottom strap for the thigh (see Fig. 9).
- 4) Engage the rear, top strap for the thigh on the pad button with the strap guide loop and pull it to the back.
 - The strap guide loop engages audibly.
- 5) Tighten the rear, top strap for the thigh (see Fig. 10).
- 6) Engage the rear, bottom strap for the lower leg on the pad button with the strap guide loop and pull it to the back.
 - The strap guide loop engages audibly.
- 7) Tighten the rear, bottom strap for the lower leg (see Fig. 11).
- 8) Push the orthosis towards the thigh and fasten the front strap for the thigh (see Fig. 12).
- 9) Fasten the front strap for the lower leg (see Fig. 13).
- 10) Tighten the orthosis straps so that the orthosis fits firmly but is comfortable to wear.
- 11) After **15 min** of active movement (e.g. slight knee bends): slightly retighten the straps in the specified sequence (see Fig. 2, straps 1-6).

5.2.5 Conversion to CI Version

The orthosis is in the delivery condition (ACL or CI version). No adaptation is required.

5.2.6 Conversion to PCL Version

- > The orthosis is in the delivery condition (ACL/CI version).
- 1) Switch the Y-hook-and-loop of the front strap for the thigh (strap 2) with the Y-hook-and-loop of the rear, bottom strap for the thigh (strap 5).
- 2) Position the orthosis correctly (see the section "Applying the ACL Version", see Fig. 7).
- 3) Fasten the straps in the modified sequence (see Fig. 3, straps 1-6).

5.3 Delivery

NOTICE

Incorrect setting

Damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ The orthosis may only be adjusted by qualified personnel.
- ▶ Do not make any improper changes to the settings.
- ▶ **Inform the patient.**

Check for proper fit of the orthosis on the leg at the time of delivery. Ensure that the centre of the patella and the centre of the joints are approximately lined up with each other.

5.4 Cleaning

INFORMATION

- ▶ Ensure proper alignment when inserting the tibia pad. The wide side of the pad has to face up.

- 1) Remove the straps and pads.
- 2) Hand-wash the straps and pads in warm water at 30 °C (86 °F) using neutral detergent.
- 3) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).
- 4) If the orthosis comes into contact with salt water or dirt, rinse with clear water and allow to air dry.

5.5 Information on Re-Use (Applies for Germany Only)

The product is approved for re-use **exclusively** in Germany.

Products that are being re-used are subject to increased strain. Re-use is only permitted up to a maximum total period of use of **2** years.

Reconditioning

Hygienic reconditioning of the product (according to the Robert-Koch Institute guidelines) is required prior to re-use.

All components subject to pronounced wear or that come into direct contact with skin are not suitable for re-use. These components are contained in the 29K110N Replacement Kit.

All components suitable for re-use must be checked for proper functionality and damage prior to re-use.

5.5.1 Replacing the Straps and Pads

Straps:

- 1) The upper and lower straps can be released from the frame using the quick-release locks.
- 2) For the straps closest to the knee, remove the strap attachment caps on the outside of the frame.
- 3) Release the plastic pin from the frame, insert the new strap and secure with the cap.

Strap and frame pads:

- The strap and frame pads are attached with hook-and-loop closures and are easy to replace.

6 Legal Information

6.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

6.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH