Otto Bock®

QUALITY FOR LIFE

50S9 Malleo Direxa Stirrup

(D)	Sprunggelenkorthese	3
(GB)	Ankle Orthosis	7
(F)	Orthèse de cheville	11
(I)	Ortesi per caviglia	15
(E)	Órtesis de tobillo	19
(P)	Ortótise da articulação do tornozelo	23
(NL)	Enkelbandage	27
(S)	Fotledsortos	31
(DK)	Ankelledsortose	35
(N)	Fotleddkortese	39
(FIN)	Nilkanivelortoosi	43
(PL)	Orteza stawu skokowego	47
(H)	bokaizületi ortézis	51
(CZ)	Ortéma kotníku	55
(RO)	Orteză pentru gleznă	59
(HR)	Ortoza skočnog zgloba	63
(BiH)	Ortoza skočnog zgloba	67
(TR)	Ayak eklemleri kemiği ortezi	71
(GR)	Νάρθηκας ποδοκνημικής	75
(RUS)	Голенистопный ортез	79
(J)	足関節装具 マレオ ディレクサ スターアップ	84

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⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION

Please read these instructions for use carefully and pay special attention to the safety information!

1 Description

1.1 *Intended Use*

The 50S9 Malleo Direxa Stirrup ankle orthosis is to be used **exclusively** for the orthotic fitting of the lower extremities.

1.2 *Indications and Effects*

The 50S9 Malleo Direxa Stirrup ankle orthosis is used for functional treatment of ligament instabilities, post-operatively after ligament ruptures, after ankle distortions, as well as for functional-conservative therapy of lateral collateral ligament ruptures. Specific indications must be determined by the physician.

The orthosis stabilises the ankle joint. Eversion and inversion of the foot joint are restricted without negatively affecting plantar and dorsal flexion of the foot.

1.3 *Contraindications*

1.3.1 **Absolute Contraindications**

None.

1.3.2 **Relative Contraindications**

In case of the following indications, consultation with a physician is required: Skin diseases/injuries, inflammatory symptoms, prominent scars on body parts to which the orthosis is to be applied that are swollen, reddened and feel excessively warm; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the orthosis is to be applied; sensory and circulatory disorders.

1.4 Safety Instructions

INFORMATION

The patient is to be instructed in the proper use/care of the product.

The initial fitting and application of the product must be carried out by trained, qualified personnel only. The daily duration of use and period of application are dependent on medical indication by the physician.

⚠ CAUTION

Risk of injury as a result of improper use. The product is designed for use on one patient only. For hygienic reasons, stockings are always to be worn under the orthosis.

An orthosis/support applied too tightly to the body can cause local pressure and, in some cases, even restrict adjacent vessels or nerves. Do not apply the product too tightly. Consult a physician immediately if you experience unusual changes (such as increase in pain).

Improper changes to the product are not permitted.

NOTICE

Risk of damage due to incorrect environmental conditions. This product is not flame-resistant. Keep the product away from flames or other heat sources.

The product should not come into contact with grease or acidic agents, unguents and lotions, as this can reduce the useful life – especially of pads.

2 Handling

2.1 Size Selection

Since the orthosis is universally applicable, no size needs to be selected.

2.2 Adapting and Applying the Product

The initial fitting and application of the orthosis must be carried out by qualified personnel.

The orthosis is delivered ready for trial fitting.

1. To apply the orthosis, completely open it and position the foot in the middle of the orthosis (Fig. 1).
2. Position the side elements medially and laterally on the ankle joints and close the lower hook and loop closure (Fig. 2).
3. Then close the upper hook and loop closure (Fig. 3).
4. If necessary, the cooling elements can be exchanged with the textile elements (Fig. 4).
Observe the corresponding instructions for use in Section 2.3.

2.3 Instructions for Use and Care

CAUTION

Risk of injury due to incorrect use of the cooling elements. The cooling elements should not come into direct contact with the naked skin, since this could result in local frost injuries/blisters. When using the cooling elements, textile material (e. g. a towel) must always be placed between the cooling element and skin.

NOTICE

Damage caused by wear and tear. Before each use, the orthosis must be checked for functional reliability and for possible wear or damage. If a closure or any other part of the orthosis shows any signs of wear, the orthosis should no longer be used.

Improper changes to the product are not permitted.

The orthosis must be used in strict conformity with the indication. The physician's instructions must be strictly observed.

Material: Side elements: polyethylene; textile pads: velour/PU; inner side of the textile: PU; hook and loop closures: nylon straps

Cooling elements: During the acute phase, patients can use the cooling elements for cooling as long as this is comfortable for them. Otherwise, the textile elements are to be used. For "(re)charging", remove the cooling elements and put them into the freezing compartment of a refrigerator for at least 6 hours (over night, if possible, at -18 °C maximum).

Cleaning: Use a slightly moistened cloth to clean the outside of the orthosis. Hand-wash the textile elements in 30 °C warm water with a standard mild detergent. We recommend closing the hook and loop closures before washing to keep them functional for a long time. Rinse thoroughly so that all residual detergent is removed. Let air dry. Do not expose to the sun and do not use oven or radiator heat for drying.

Disposal after usage: The material used in the product is recyclable and should be disposed of accordingly.

3 Further Usage Restrictions

The product is designed for use on one patient only.

The daily duration of use and period of application are dependent on the medical indication.

4 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

5 CE Conformity

This product meets the requirements of the 93/42/EEC guidelines for medical devices. This product has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH