











C-Leg Protector 4X860

DE Gebrauchsanweisung (Fachpersonal)	4
EN Instructions for use (qualified personnel)	11
FR Instructions d'utilisation (Personnel spécialisé)	17
IT Istruzioni per l'uso (Personale tecnico specializzato)	25
ES Instrucciones de uso (Personal técnico especializado)	32
PT Manual de utilização (Pessoal técnico)	39
NL Gebruiksaanwijzing (Vakmensen)	46
SV Bruksanvisning (Fackpersonal)	53
DA Brugsanvisning (Faguddannet personale)	59
NO Bruksanvisning (Fagpersonell)	66
FI Käyttöohje (Ammattihenkilöstö)	73
PL Instrukcja użytkowania (Personel fachowy)	79
CS Návod k použití (Odborný personál)	87
HR Upute za uporabu (Stručno osoblje)	93
TR Kullanma talimatı (Uzman personel)	100
EL Οδηγίες χρήσης (Τεχνικό προσωπικό)	106
RU Руководство по применению (Квалифицированный персонал)	114
ZH 使用说明书 (专业人员)	121
KO 사용 설명서 (전문기술요원)	127

INFORMATION

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INFORMATION

Last update: 2015-02-05

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

The "C-Leg Protector" is referred to simply as the protector/product below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.


Instruct the patient in the proper use and care of the product. The product may not be transferred to the patient without prior instruction.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product Description

2.1 Design

The product consists of the following components:

1. Protector main unit
2. Upper closure (one side always stays on the main unit)
3. Lengthwise closure
4. Foot cuff for protector
5. Shield insert

2.2 Function

The protector provides a means of enclosing the C-Leg electronic knee joint, including the tube adapter, to enhance its appearance and function. The housing of the knee joint is protected against external influences by the protector, e. g. scratching and surface contamination. The protector is not required for the proper and safe functioning of the joint.

2.3 Combination Possibilities

Prosthetic components that are suited for combination with the product are listed below.

Knee Joint Combination Possibilities

Designation	Reference number
C-Leg	3C98-3/3C88-3

Prosthetic Foot Combination Possibilities

Prosthetic foot		Foot cuff		
Designation	Reference number	Size S	Size M	Size L
Triton	1C60	24-25	26-28	29-30
Triton Vertical Shock	1C61	24-25	26-28	29-30
Triton Low Profile	1C63	24-25	26-28	29-30
Triton Heavy Duty	1C64	24-25	26-28	29-30
Trias	1C30	23-25	26-28	29-30
C-Walk	1C40	24-25	26-28	29-30
Dynamic Motion	1D35	23-25	26-28	29-30
Adjust	1M10	23-25	26-28	29-30
Axtion	1E56	23-25	26-28	29-31
Lo Rider	1E57	24-25	26-28	29-31

INFORMATION

The foot cuff is **not** compatible with the slim footshells.

3 Application

3.1 Indications for use

The housing of the knee joint is protected against external influences by the protector, e. g. scratching and surface contamination.

3.2 Area of Application

The product is intended **exclusively** for C-Leg users.

In principal, the product can be used with a rotation axis-floor-measurement* of 430 mm - 560 mm.

An exception is fittings with the 1C61 Triton Vertical Shock prosthetic foot. Due to the height of the foot, the fitting area is restricted to 470 mm - 560 mm.

*Position of the alignment reference point of the C-Leg (rotation axis): 20 mm above the medial tibial plateau.

Rotation axis-floor-measurement measured while standing (barefoot)

INFORMATION

For smaller users that require the Protector to be substantially shortened, the foot cuff material can extend forwards into the area of the shield insert. The rotation axis-floor-measurement at which this may occur depends on the foot cuff size being used.

3.3 Conditions of use

The product was developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

3.4 Qualification

The fitting of patients with the product may only be carried out by qualified personnel (e.g. a prosthetist).

4 Safety

4.1 Explanation of Warning Symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

4.2 Structure of the safety instructions

CAUTION

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General Safety Instructions

NOTICE

Improper assembly/disassembly of the product

Damage to the product.

- ▶ Observe the process steps, instructions and required tools in the section "Preparation for use".
- ▶ Instruct the patient in the proper disassembly/ assembly of the product.

NOTICE

Independent user changes or modifications made to the product

Damage to the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.

5 Scope of Delivery

Scope of delivery for 4X860=*

- 1 pc. C-Leg Protector (without shield insert)
- 1 pc. 4P860=R lengthwise closure
- 1 pc. 4P860=U upper closure
- 1 pc. 4P880=* foot cuff for protector
- 1 pc. 647G1113 instructions for use (qualified personnel)
- 1 pc. 647G1114 instructions for use (user)

Accessories

The following components are not included in the scope of delivery and must be ordered separately:

- 1 pc. 4P863=* shield insert

6 Preparation for Use

Required Tools

Designation	Reference number
Vibrating/band saw	-
Knife	-
Loctite® 241 Thread Locking Compound	636K13
Crepe Adhesive Tape	627B6=15
Torque Wrench	-
Tape Measure	-

6.1 Adapting the Protector

6.1.1 Shortening the Protector – Preliminary Shortening

- 1) Remove the prosthesis from the residual limb.
- 2) Remove the shoe from the prosthetic foot.
- 3) Loosen the foot connection cap of the prosthetic foot.
- 4) Check whether the closures on the protector are properly fastened and the protector is not twisted.
- 5) Lay the protector next to the prosthesis.
INFORMATION: Make sure that the knee plateaus of the knee joint and protector are at the same level (see Fig. 1).
- 6) Mark the protector distally at the level of the foot contact surface (see Fig. 2).
INFORMATION: For tall patients, shortening the protector may not be required.
- 7) Make another mark 1 cm above the existing mark.
- 8) Apply crepe adhesive tape proximally all around at this mark.
- 9) Use a vibrating or band saw to shorten the protector along the bottom edge of the crepe adhesive tape.

6.1.2 Shortening the Protector – Final Shortening

- 1) Remove the crepe adhesive tape and lengthwise closure from the protector after preliminary shortening.
- 2) Open the upper closure on the protector.
- 3) Spread open the protector and apply it to the knee joint (see fig. 4).
INFORMATION: The retainers of the protector (at the side in the upper closure area) have to engage on the frame element of the knee joint (see fig. 5).
- 4) Fasten the upper closure on the protector.
→ Slightly squeeze the protector together with both hands.
- 5) Depending on the size of the chosen foot cuff, mark the following distances in at least 3 places on the protector (see fig. 6):
INFORMATION: The distances have to be measured respectively from the foot contact surface, creating a cut edge parallel to the foot contact surface. This makes for an optimum overlap of the foot cuff over the protector main unit.
 - Size S: 110 mm
 - Size M: 120 mm
 - Size L: 125 mm
- 6) Remove the protector.
- 7) Fasten the upper closure (see fig. 14) and lengthwise closure (see fig. 7).
- 8) Tape off the protector with crepe adhesive tape above the mark.

NOTICE

Loss of stability due to shortening the protector too much

No hold of the protector on the knee joint.

- At the most, the protector may be shortened to just before the marked rib (see fig. 3).

6.2 Finishing the Protector

- 1) Using a vibrating or band saw, shorten the protector along the bottom edge of the crepe adhesive tape (see Fig. 8).

INFORMATION: Make sure that the closures are not shifted while sawing.

- 2) Deburr the cut edge.

Additional Instruction for Using a Torsion Adapter

- Using suitable tools, remove the profile ribs on the interior of the protector (see Fig. 9).

6.3 Installing the Protector and Foot Cuff

6.3.1 Installing the Foot Cuff

- 1) Remove the prosthetic foot from the tube adapter.
In doing so, unscrew the two set screws that are screwed in the furthest but are not opposite each other. Mark the two set screws.
- 2) Pull the foot cuff over the tube adapter (see Fig. 10).
- 3) Install the prosthetic foot on the tube adapter (see Fig. 11).
Apply Loctite® 241 to the two set screws and screw them in.
- 4) Position the clips of the foot cuff on the footshell, aligning the seam to posterior.
- 5) Attach the clips in the footshell (see Fig. 12).

INFORMATION: When using the 1C40 C-Walk, the clips that do not engage have to be removed.

6.3.2 Protector Installation

- 1) Remove the crepe adhesive tape from the protector after final shortening.
- 2) Open or remove the closures.
- 3) Spread open the protector.
- 4) Thread the protector into the foot cuff from anterior (see fig. 13).
- 5) Apply the protector to the knee.

INFORMATION: The retainers of the protector have to engage on the frame element.

- 6) Fasten the upper closure (see fig. 14) and lengthwise closure (see fig. 7) on the protector.

NOTICE

Loss of stability due to inadequately installed retainers

No hold of the protector on the knee joint.

- Ensure that all retainers and closures are properly engaged and/or fastened.

6.4 Installing and removing the shield insert

- Slide the shield insert into the recess of the applied protector from anterior and let it engage (see fig. 15).

NOTICE

Loss of stability due to inadequately installed protector main unit

No hold of the shield insert on the protector.

- Ensure that the protector has been installed properly and that all closures are fastened.
- By opening the closures and removing the protector from the knee joint, the shield insert loses its hold and can be removed.

6.5 Functional Test

After assembly, the following points have to be checked:

- Does the foot cuff overlap the protector by approx. 2 cm?
- Is the function of the rotation adapter impaired?
Material may be removed from the protector to improve the functionality.

INFORMATION: If material is removed, the stability/flexibility of the protector on the knee joint may be impaired.

INFORMATION

Charging the C-Leg batteries with the protector installed

It may be difficult for the prosthesis user to open the charging receptacle cover through the recess in the lengthwise closure with his/her fingers. In this case, the battery charger plug may be used to pry open the cover.

7 Maintenance

- During the standard inspection of the prosthesis components, also inspect the product.

7.1 Cleaning and Care

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- Only clean the product with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).
- 1) Clean the product with a damp, soft cloth.
 - 2) Dry the product with a soft cloth.
 - 3) Allow to air dry in order to remove residual moisture.

8 Legal Information

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as class I accessories according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

8.3 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

9 Appendices

9.1 Symbols Used

9.1.1 Symbols on the product

 Legal manufacturer

 Declaration of conformity according to the applicable European directives

 Lot number

9.1.2 Symbols on the shield insert

 Legal manufacturer

 Declaration of conformity according to the applicable European directives

 Lot number

9.2 Technical data

Ambient Conditions	
Storage and transport in original packaging	-25 °C/-13 °F to 70 °C/158 °F
Storage and transport without packaging	-10 °C/14 °F to 60 °C/140 °F max. 100% relative humidity, non-condensing
Operation	-10 °C/14 °F to 60 °C/140 °F max. 100% relative humidity, non-condensing
General information	
Reference number	4X860=* + 4P863=*
Product service life	Wear part subject to normal wear and tear
Protector weight (with closures)	450 g
Foot cuff weight	60 g
Shield insert weight	63 g

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This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Wien · Austria
T +43-1 523 37 86 · F +43-1 523 22 64
info.austria@ottobock.com · www.ottobock.com