Otto Bock®

QUALITY FOR LIFE


CE

50P10 Manu Immobil / 50P11 Manu Immobil Long

(D)	Handlagerungssorthesen	3
(GB)	Wrist positioning Orthoses	5
(F)	Orthèses d'immobilisation de la main	8
(I)	Ortesi per mano e polso	10
(E)	Muñequera	13
(P)	Ortóteses de posicionamento da mão	16
(NL)	Handorthesen	18
(S)	Handstödsortoser	21
(DK)	Håndortoser	23
(PL)	Orteza stabilizująca dłoń i przedramię	26
(H)	Kartámasztó ortézis	28
(CZ)	Imobilizační ortézy zápěstí	31
(TR)	El Yataklama Ort ezleri	33
(GR)	Πηχεοκαρπικοί νάρθηκες	36
(RUS)	Иммобилизирующие ортезы кисти	38
(中国)	手臂固定式/手臂固定式加长托架矫形器	41

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1 Intended Use

English

The Manu Immobil 50P10 and Manu Immobil Long 50P11 Wrist positioning Orthoses are to be **exclusively** used for the orthotic fitting of the hand or wrist in accordance with the indications. They can be used for post-operative and post-traumatic immobilization as well as for occupational therapy.

The Manu Immobil Long 50P11 wrist positioning orthosis offers additional immobilisation of the finger joints.

2 Indications and Effects

The orthoses provide for a high degree of immobilization of the hand, wrist (applies to 50 P10 and 50P11) and finger joints (only applies to 50P11) after fractures and luxations or following the wrist, carpal and metacarpal operations. For treatment of acute irritations (such as carpal tunnel syndrome), wrist arthroses or severe distortions. Indication must be determined by the physician. The hand or wrist and finger joints (only applies to 50 P11) are highly immobilized in a functional position.

3 Selection of Orthosis Sizes

The **50P10 Manu Immobil Wrist Orthosis** is available in one universal size. It fits right and left sides.

The **50P11 Manu Immobil Long Wrist Orthosis** is available in two sizes (M and L) and fits right and left sides.

4 Material and Construction

The wrist orthoses consist of polypropylene with an integrated steel joint bar as well as attached pad of SpaceTex (polyester, polyamide) and tie-down straps (polyamide) with hook and loop closure.

5 Trial Fitting, Application and Delivery of the Orthosis

The initial fitting and application of the Manu Immobil and Manu Immobil Long Wrist positioning Orthoses must only be carried out by qualified personnel. The orthoses are delivered closed and ready for trial fitting and are already preset in a functional position.

If necessary, the orthopedic technician can perform individual adjustment ($\pm 20^\circ$) in the area of the wrist and metacarpophalangeal joints. The stamped lines on the underside of the orthosis mark the area where the orthosis can be cold-shaped.

Note:

Use a level support surface for adjusting the orthosis. Shape carefully, step by step, applying the pressure directly on the inner, centrally positioned steel splint. Repeated adjustment should be avoided.

To apply the orthosis, adjust the fixation straps with adjustable direction of closure to the largest extent possible. Slide the hand into the wrist orthosis from the back until the hand rests completely on the hand pad. Subsequently pull on the fixation straps in the following order and close the hook and loop closures (fig. 2 + 3):

- 1) Wrist strap
- 2) Back-of-the-hand strap
- 3) Finger strap
- 4) Forearm strap

The fixation straps, especially the back-of-the-hand strap should be applied firmly, without constriction.

Direction of strap closures

The Wrist positioning Orthoses are delivered preset for right-hand fitting. If need be, the direction of strap closures can be turned (from inside to outside or from outside to inside) as follows. First open all hook and loop closures. Remove the pad inclusive the three Velcro points (fig. 4). Remove, turn them and re-assemble the straps. Make sure that the plastic loops are located on the underside of the orthosis.

Length of straps

All straps can be trimmed if necessary. For this purpose undo the hook and loop closures of the plastic loop of the respective strap, trim the strap and re-attach the plastic loop.

Positioning of the finger strap (with 50P11 only)

If necessary, individually position the finger strap through a wide slot on the underside of the orthosis. Remove the pad of the orthosis, move the finger strap to the desired position and re-attach the pad.

Upon delivery of the orthosis to the patient, the orthosis must be tested for correct fit. Upon delivery of the product, the patient is to be instructed in the operation of the 50P10 Manu Immobil and 50P11 Manu Immobil Long wrist positioning orthoses.

6 Side Effects

Attention!

All appliances, such as supports and orthoses, applied too tightly to the body, can cause local pressure and, in rare cases, restrict blood circulation or nerves.

7 Contraindications

Attention!

In case of the following clinical pictures, the application and use of this type of aid is indicated only after consultation with the physician.

- 1) Skin diseases and injuries of body parts to be fitted, especially in case of inflammatory symptoms. This also applies to prominent scars that are swollen, reddened and feel excessively warm when touched.
- 2) Lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the appliance will be applied.
- 3) Sensory disorders in the arm/hand.

8 Instructions for Use and Care

The orthosis must be used in strict conformity with the indication. The patient's acceptance and cooperation play an important part. It is essential to instruct the patient in the proper use of the orthosis! The physician's instructions must be strictly observed.

The daily duration of use and total period of application must be determined by the attending physician.

Patients should be instructed to consult a physician immediately if they experience unusual changes (such as an increase in pain).

Attention!

The orthosis is not flame-resistant. Patients should be instructed to keep the product away from flames or other heat sources.

Before each use, the orthosis must be checked for functional reliability and for possible wear or damage. If a hook and loop closure or any other part of the orthosis show any signs of wear, the orthosis should no longer be used.

Attention!

The patient should be instructed not to wear the orthosis during nuclear spin tomography examination.

The product should not come into contact with grease or acidic agents, unguents and lotions.

Improper changes to the product are not permitted.

You can remove the pad and Velcro straps and wash the orthosis by hand at 30°C using a standard mild detergent. We recommend closing the hook & loop closures before washing to prevent any damage to them. To avoid skin irritations, rinse the orthosis thoroughly so that all soap residue is removed. Do not expose to the sun and do not use oven or radiator heat for drying. Let the support air dry.

9 Restrictions for Use

The 50P10 Manu Immobil and 50P11 Manu Immobil Long Wrist positioning orthoses are designed to be used on **one** patient. The daily duration of use and period of application are dependent on medical indication. Wounds are to be protected with a dressing.

Attention!

Parts to be fitted and those parts that come directly into contact with the skin can cause functional and hygienic risks if the orthosis is used by another person.

10 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

11 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH