



28K20, 28K21 Genu Arthro


DE Gebrauchsanweisung	7
EN Instructions for Use	12
FR Instructions d'utilisation	17
IT Istruzioni per l'uso	23
ES Instrucciones de uso	28
PT Manual de utilização	33
SV Bruksanvisning	39
PL Instrukcja użytkowania	44
TR Kullanma talimatı	49
EL Οδηγίες χρήσης	54

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INFORMATION

Last update: 2013-11-15

- Please read this document carefully.
► Follow the safety instructions.

These instructions for use provide important information on the fitting and application of the Genu Arthro knee joint orthosis, comprising the 28K20 Transfemoral Module and the 28K21 Transtibial Module.

2 Product Description

2.1 Components

The orthosis comprises two main modules with the following components:

28K20 Transfemoral Module (see Fig. 1)

No.	Component	No.	Component
1	Tube frame with upper joint section and end cap	3	Thigh pad
2	Belt guiding element for thigh	4	Thigh straps and belt pads

28K21 Transtibial Module (see Fig. 2)

No.	Component	No.	Component
1	Tube frame with lower joint section, slide washers, joint screws and end cap	6	Knee support with 3 different pressure pads
2	Belt guiding element for calf	7	Bearing washer

No.	Component	No.	Component
3	Calf pads	not illustrated	Key for assembling joint
4	Calf straps and belt pads	not illustrated	Torx key for adjusting tube frame
5	Foot section for anti-migration stability		

The bearing washer (no. 7) is available in 4 versions distinguishable by their colour.

Thickness [mm]-0.05	0.4	0.5	0.6	0.7
Colour	blue	natural	black	red

2.2 Design

The orthosis is available in four different sizes. It is possible to provide a wide spectrum of patients with an individual fitting thanks to the ability to combine any size of the transfemoral module and the transtibial module.

The required sizes can be determined using the size table.

3 Application

3.1 Intended Use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

3.2 Indications

- Single-compartment osteoarthritis of the knee

Indications must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

Not known.

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the medical device will be applied; sensory and circulatory disorders in the leg area.

3.4 Effects

The orthosis provides relief for the affected compartment by means of the 3-point principle and stabilises the knee joint.

4 Safety

4.1 Explanation of Warning Symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.

4.2 General Safety Instructions

CAUTION

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

CAUTION

Contact with heat, embers or fire

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.
- ▶ **Inform the patient.**

CAUTION

Reuse on other patients and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

NOTICE

Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.
- ▶ **Inform the patient.**

5 Preparation for Use

5.1 Size Selection

- 1) Determine the size of the transfemoral module (28K20) from the distance between the crotch height and medial tibial plateau (see size table – measurement A).
- 2) Determine the size of the transtibial module (28K21) from the distance between the medial tibial plateau and 5 cm above the lateral ankle (see size table – measurement B).

5.2 Assembly

CAUTION

Incorrect alignment, assembly or adjustment

Injuries may result from incorrectly assembled, incorrectly adjusted or damaged components

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.
- ▶ **Inform the patient.**

> **Required tools and materials:** key for assembling joint, Torx key for adjusting tube frame, Loctite® (636K13), torque wrench.

- 1) Remove the joint screw from the lower joint section using the key provided.
- 2) Remove the upper slide washer (see Fig. 3).

- 3) Insert the upper joint section (see Fig. 4).
- 4) Position the slide washer that you have removed so that it compensates sideways play in the knee joint of the orthosis.
- 5) Insert the bearing washer (2) between the knee joint (1) and joint nut (3) of the transfemoral and transtibial module (see Fig. 5).

INFORMATION: The thickness of the bearing washer that needs to be inserted depends on the amount of play.

- 6) Reinstall the joint screw and secure it with Loctite® if applicable.
- 7) Loosen the screws in the end cap of the transtibial module.
- 8) Insert the foot section into its channel (see Fig. 6).

INFORMATION: Use the distance between the patient's knee and the floor as a reference for the initial positioning of the foot section.

- 9) Retighten the screws in the end cap of the transtibial module (**tightening torque 3 Nm**).
- 10) Attach the straps with ergonomic closures to the dorsal side of the orthosis.

6 Settings

6.1 Adjusting the Orthosis

CAUTION

Improper moulding or application

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.
- ▶ **Inform the patient.**

- 1) Determine the knee joint's compromise pivot point according to Nietert.
- 2) Position the orthosis joint in accordance with this (see Fig. 7 and see Fig. 8).
- 3) **If necessary:** Readjust the pre-settings of the foot section.
- 4) Adapt the strap guiding elements to the contours of the thigh (ventrally) and the calf (dorsally).
- 5) Slightly loosen the screws in the tube clamp at the end caps and at the joint (see Fig. 9), in order to bring the frame to as neutral a position as possible (free of pressure) next to the leg.
- 6) Adapt the frame to the contour of the leg by pulling or pressing lightly (see Fig. 10 and see Fig. 11).
- 7) **If the adjustments made are insufficient:** Adjust the tube frame by hand, e.g. using a ledge.

NOTICE! Functional risks are associated with, for example, repeated cold forming and the resulting reduction in breaking elongation (strain hardening).

- 8) Retighten the screws (**tightening torque 3 Nm**).
- 9) Hook on the lower quick-release fastener on the thigh ventrally.
- 10) Hook on the upper quick-release fastener on the calf ventrally.
- 11) Hook on the lower quick-release fastener on the calf ventrally.
- 12) Hook on the upper quick-release fastener on the thigh ventrally.
- 13) **If necessary:** Remove the Y-shaped hook-and-loop strap from the belt and cut down the belt until there is approx. 5 cm excess material to the quick-release fastener (see Fig. 13).
- 14) Put the Y-shaped hook-and-loop strap back on and fasten the belt in the same order.
- 15) Cut the belt pads accordingly (see Fig. 14) and insert them.

6.2 Adjusting the Corrective Effect

- 1) Insert the support with the first, thin pad into the joint in order to create the axis-correcting pressure.

2) **If necessary:** Insert the second or third pad.

Another possibility to build up corrective pressure is to use the tube frame.

- 1) Loosen the thigh straps.
- 2) Loosen the screws in the tube clamp at the end caps and at the joint.
- 3) Bring the end caps to within approx. 1-2 cm from the leg by pulling lightly.
- 4) Retighten the screws in the tube clamp (**tightening torque 3 Nm**).
- 5) Fasten the straps.
- 6) **If necessary:** Check the corrective effect using the LASAR Posture.

7 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

7.1 Application

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.
- ▶ **Inform the patient.**

CAUTION

Using orthosis without padding

Risk of local pressure points, superficial injuries or skin irritation due to skin contact with rigid or sharp orthosis components

- ▶ Only use the orthosis with undamaged padding.
- ▶ **Inform the patient.**

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear and damage prior to each use.
- ▶ Advise the patient that the product must not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

- > Open the quick-release fasteners on the orthosis.
- > The patient is sitting on the front edge of a chair.
- 1) Place the orthosis on the leg. Align the orthosis joint.
- 2) Position the foot section in the shoe (see Fig. 15).
- 3) Hook on the lower quick-release fastener on the thigh ventrally.
- 4) Hook on the upper quick-release fastener on the calf ventrally.

- 5) Hook on the lower quick-release fastener on the calf ventrally.
- 6) Hook on the upper quick-release fastener on the thigh ventrally.

7.2 Cleaning

Clean the orthosis regularly:

- 1) Remove the straps and pads from the orthosis.
- 2) Close all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at 30 °C (86 °F) using neutral detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

7.3 Care

► **If necessary:** Lubricate the orthosis joint with a dry lubricant (e.g. 633T2 Teflon® Spray).

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal Information

9.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

9.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH