1C10 Terion


DE	Gebrauchsanweisung (Fachpersonal)	5
EN	Instructions for use (qualified personnel)	12
FR	Instructions d'utilisation (Personnel spécialisé)	19
IT	Istruzioni per l'uso (Personale tecnico specializzato)	27
ES	Instrucciones de uso (Personal técnico especializado)	35
PT	Manual de utilização (Pessoal técnico)	43
NL	Gebruiksaanwijzing (Vakmensen)	50
SV	Bruksanvisning (Fackpersonal)	58
DA	Brugsanvisning (Faguddannet personale)	65
NO	Bruksanvisning (Fagpersonell)	72
FI	Käyttöohje (Ammattihenkilöstö)	79
PL	Instrukcja użytkowania (Personel fachowy)	86
HU	Használati utasítás (szakszemélyzet)	94
CS	Návod k použití (Odborný personál)	101
RO	Instrucţiuni de utilizare (Personal de specialitate)	108
HR	Upute za uporabu (Stručno osoblje)	116
SL	Navodila za uporabo (Strokovno osebje)	124
SK	Návod na používanie (Odborný personál)	131
BG	Инструкция за употреба (Специалисти)	138
TR	Kullanma talimatı (Uzman personel)	146
EL	Οδηγίες χρήσης (Τεχνικό προσωπικό)	153
RU	Руководство по применению (Квалифицированный персонал) ...	161
JA	取扱説明書 (有資格担当者)	169
ZH	使用说明书 (专业人员)	176
KO	사용 설명서 (전문기술요원)	182

1

2

7

INFORMATION

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INFORMATION

Last update: 2015-08-31

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The 1C10 Terion prosthetic foot is suitable for walking on various surfaces and for a large range of speeds.

The functional properties of the prosthetic foot are achieved through the carbon spring element and a functional foam base. The prosthetic foot permits perceptible plantar flexion at heel strike and a natural rollover movement. Stored energy is released by the spring element.

1.2 Combination Possibilities

This prosthetic component is part of the Ottobock modular system. As a result of its design, it may be combined with other prosthetic components of the

modular system. You can find possible restrictions in this section. Contact the manufacturer with questions.

2 Application

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of Application

The product is recommended for mobility grade 2 (restricted outdoor walker) and mobility grade 3 (unrestricted outdoor walker).

Size [cm]	22 to 23	24 to 25	26 to 28
Max. body weight [kg]	75	100	125

2.3 Environmental Conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +60 °C
Relative humidity: no restrictions
Chemicals/moisture: fresh water, salt water, chlorine water, soapsuds, sweat, urine
Solids: dust, sand, highly hygroscopic particles (e. g. talcum)
Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage (see Page 18).

2.4 Service Life

This product has been tested by the manufacturer with 2 million load cycles according to ISO 22675. Depending on the patient's activity level, this corresponds to a service life of 2 to 3 years.

3 Safety

3.1 Explanation of Warning Symbols

CAUTION	Warning regarding possible risks of accident or injury.
NOTICE	Warning regarding possible technical damage.

3.2 General Safety Instructions

CAUTION

Exceeding the service life and reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded (see section "Service life").
- ▶ Only use the product for a single patient.

CAUTION

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the prosthetic components according to the MOBIS classification (see section "Area of application").

CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved according to the section "Combination possibilities".
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

NOTICE

Mechanical damage to the product

Change in or loss of functionality due to damage

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

NOTICE

Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions (see the section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Signs of changes in or loss of functionality during use

Reduced spring effect or decreased forefoot resistance are indications of loss of functionality.

4 Scope of Delivery

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

1C10 Terion				
Fig.	Item	Quantity	Designation	Reference number
–	–	1 (■)	Instructions for Use	647G972
–	–	1 (■)	Prosthetic Foot	–
–	–	1 (■)	Footshell	2C12=*
–	–	1 (■)	Connection cap	2C13=*
–	–	1 (■)	Toe insert	2Z661=*
–	–	1 (■)	Protective Sock	SL=Spectra-Sock

5 Preparation for Use

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

NOTICE

Grinding the prosthetic foot or footshell

Premature wear resulting from damage to the product

- Do not grind the prosthetic foot or footshell.

5.1 Alignment

INFORMATION

There is a plastic adapter cover on the proximal connection of the prosthetic foot. This protects the connecting section from scratches during the alignment and trial fitting.

- Remove the adapter cover before the patient leaves the workshop/fitting area.

5.1.1 Applying/removing the footshell

Inserting the foot into the footshell

- > **Recommended tools:** 2C100 changing device for footshell
- 1) Pull the Spectra sock over the prosthetic foot, avoiding any wrinkles.
- 2) Slide the prosthetic foot into the footshell. Insert the carbon spring into the slot on the toe insert (see fig. 4).
- 3) Push the prosthetic foot into the heel area of the footshell until it locks into place (see fig. 6).
- 4) Place the connection cap on the footshell (see fig. 7).

Removing the footshell

- > **Recommended tools:** 2C100 changing device for footshell
- 1) Push the heel area of the footshell to the rear and pull the prosthetic foot up.
- 2) Remove the prosthetic foot from the footshell.

5.1.2 Bench Alignment

TT bench alignment

Bench alignment process	
	Recommended tools: 662M4 goniometer, 743S12 heel height measuring apparatus, 743A80 50:50 gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. assembly or 743A200 PROS.A. assembly)
Perform the alignment and assembly of the prosthetic components in the alignment apparatus according to the following specifications:	
①	Heel height: Effective heel height (shoe heel height – sole thickness of fore-foot) + 5 mm

Bench alignment process	
2	Frontal alignment: Alignment reference line on a–p markings of the connection cap/footshell (see fig. 2)
3	Sagittal alignment: Alignment reference line on m–l markings of the connection cap/footshell (see fig. 2)
4	Connect the prosthetic socket and prosthetic foot using the chosen adapters. Be sure to follow the adapter instructions for use during adjustment and installation.
5	Sagittal alignment: Determine the centre of the prosthetic socket with the 50:50 gauge. Align the prosthetic socket centrally to the alignment reference line. Socket flexion: Residual limb flexion of the patient + 5°
6	Observe the abduction or adduction position.

TF bench alignment

- Observe the information in the prosthetic knee joint instructions for use.

5.1.3 Static Alignment

- Use the L.A.S.A.R. Posture for static alignment optimisation of the prosthesis, if available.
- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) may be requested from Ottobock.

5.1.4 Dynamic Trial Fitting

- Adapt the setup of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) so that correct heel contact and optimal roll-over are ensured.
- **TT fittings:** Make sure that physiological knee movement is achieved after the heel makes impact and the leg begins to bear weight.
- Remove the adapter cover from the prosthetic foot after completion of the dynamic trial fitting and gait training exercises.

5.2 Optional: Installing the foam cover

The connection cap serves as a detachable connection between the cosmetic foam cover and footshell, and snaps precisely into the proximal brim of the footshell.

- > **Required materials and tools:** 634A58 Isopropyl Alcohol, 636N9 Contact Adhesive or 636W17 Plastic Adhesive
- 1) Cut the foam cover to length according to the patient's measurements, allowing for distal compression of about **1 cm** for transtibial and of about **4 cm** for transfemoral covers.
 - 2) Pull the foam cover over the prosthesis.
 - 3) Join the connection cap to the footshell and mount the prosthetic foot.
 - 4) Mark the outer contours of the connection cap on the distal face of the foam cover blank.
 - 5) Dismount the prosthetic foot and remove the connection cap from the footshell.
 - 6) Clean the connection cap with isopropyl alcohol.
 - 7) Glue the connection cap to the face of the foam cover as marked and allow the glued joint to set (approx. **10 minutes**).
 - 8) Mount the prosthetic foot and adapt the exterior cosmetic cover shape. Take into account compression due to cosmetic stockings or SuperSkin.

6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

9 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

10 Technical data

Sizes [cm]	22	23	24	25	26	27	28
Heel height [mm]	5 ±5						
Footshell	slim		normal				
System height with foot-shell [mm]	43				49		52
Structural height [mm]	61				67		
Product weight without footshell [g]	245		302		343		
Product weight with foot-shell [g]	411	439	506	528	598	602	661
Max. body weight [kg]	75		100		125		
Mobility grade	2 + 3						

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Otto Bock HealthCare GmbH