

## Neptun Super Soft Height Adjustable Shower Stool



This version of our famous Neptun Shower Stool is fitted with a very soft PU foam seat and backrest (Backrest is an accessory)

Made from stainless steel  
(with telescopic legs made of aluminium)

Very easy to clean

Easily height adjustable (43-64 cm)

The backrest can be mounted or removed without the use of tools.

The backrest is depth adjustable in two steps

The stool is available in two widths (41 and 50 cm)

Seat and back can be cleaned in disinfecting bedpan washers



⇐ Soft PU seat that snaps on easily onto the frame.

The seat has grooves and holes for draining excess water from the seat.



Extra soft seat for even higher comfort.



⇐ Seat height is adjusted with clips.

The backrest⇒ Is mounted with clips, and adjustable in two different depths. The foam piece on the back can be removed for cleaning.



**HMN** since 1934

**Technical data:****Neptun Super Soft, Height Adjustable Shower stool, with extra soft seat**

	Shower Stool 41cm	Shower Stool 50 cm	Backrest
Seat height	43-64 cm	43-64 cm	-
Total height (upper armrest)	56-77 cm	56-77 cm	Upper edge of backrest 81-102 cm
Width between Armrests	41 cm	50 cm	-
Size Foam Seat	B40xD36,5xH4 cm	B40xD36,5xH4 cm	-
Width at the floor, approximately		61-65 cm	-
Depth at the floor, approximately	43-49 cm	43-49 cm	-
Measurements, backpiece – PU part	-	-	Width lower edge 34 cm Width upper edge 24 cm Height 27 cm
Adjustable Depth – Backrest	-	-	2 cm from edge of seat 4 cm from edge of seat
Weight	Approx 5 kg	Approx 5,9 kg	Approx 1 kg
Material: Seat /backrest	PU Foam	PU Foam	PU Foam
Maintenance	Normal cleaning agents at max. 80 °C	Normal cleaning agents at max. 80 °C	Same agents as the bench itself
Stackable	Yes	Yes	-
Maximum userwheight	150 kg	150 kg	150 kg
Frame	Stainless steel, powder coated, white	Stainless steel, powder coated, white	Stainless steel, powder coated, white
Teleskope legs	Aluminium (stainless) height adjustable, 8 levels	Aluminium (stainless) height adjustable, 8 levels	-
Rubber feet for legs	Grey plastic, with draining hole	Grey plastic, with draining hole	-
Description:		Art. no	
Neptun Super Soft height adjustable shower stool 41 cm		312002	
Neptun Super Soft height adjustable shower stool 50 cm		312003	
Neptun Super Soft Backrest		312021	
Spare parts:			
Clips for height adjustment and mounting of backrest		310065	
PU foam Seat for Neptun Super Soft		800300	
Teleskopic leg (1 piece with clip, rubber foot and bush)		312033	
Rubber foot		310017	
Bush		310021	
PU foam backrest piece Neptun Super Soft		800301	



This product is CE marking which means the product complies with EU legislations and enables the free movement of the product in the European market.

07.10.15



**HMN** since 1934