

Neptun Shower Stool, height adjustable, ALUMINUM



Frame and legs produced in stainless aluminum.

Comfortable seat and backrest.

Very easy to clean.

Easy to adjust height (43-64 cm)

Backrest (=accessory) can be put on and taken off without using tools.

The backrest can be adjusted in two depths.

The shower stool comes in two different widths (41 and 50 cm)



⇐ Air filled plastic seat, which you can click off (for cleaning a.s.o.). The seat is with a nobly surface and with water grooves.

⇒ The frame is designed to make it very easy to clean.



⇒ The backrest is put on with two clips. It can be adjusted in two depths. The air filled plastic part can be clicked off (for cleaning a.s.o.)

⇐ The Height is adjusted with clips



Technical data:
Neptun Shower Stool, Height adjustable, ALUMINUM

	Shower Stool 41 cm	Shower Stool 50 cm	Backrest
Seat heights	43-64 cm	43-64 cm	-
Total heights (at the armrest)	56-77 cm	56-77 cm	At the top of the backrest 79-100 cm
Width between armrests	41 cm	50 cm	-
Seat size.(size of the plastic seat)	W41xD36,5xH4,5 cm	W41xD36,5xH4,5 cm	-
App. Widths at floor	54,5-58,5 cm	60,5-64,5 cm	-
App. Depths at floor	44-52 cm	44-52 cm	-
Size of the plastic part of the backrest	-	-	Width bottom 33 cm Width top 27 cm Height 20 cm
Depth backrest (can be adjusted in two different depths)	-	-	2 cm from seat edge and 4 cm from seat edge
Weight	3,4 kg	3,5 kg	1,3 kg
Seat/back plastic part material	Polythene plastic	Polythene plastic	Polythene p.
Cleaning	Ordinary cleaning detergents at max. 80 degrees.	Ordinary cleaning detergents at 80 degrees.	Same as the shower stool
Stackable	Yes	Yes	-
Max. User load	130 kg	130 kg	150 kg
Frame	Aluminum. Powder lacquered	Aluminum. Powder lacquered	Stainless Steel Powder Lacquered
Telescopic legs	Aluminium. Height adjustable in 8 positions.	Aluminium. Height adjustable in 8 positions.	-
Rubber feet	Grey plastic. Hole in bottom for outflow of water.	Grey plastic. Hole in bottom for outflow of water	-

Description:	Art. no
Neptun Height adjustable shower stool, ALUMINUM width 41 cm	312006
Neptun Height adjustable shower stool, ALUMINUM width 50 cm	312007
Neptun backrest	312020
Spare parts:	
Clips for height adjustment (the same one for backrest)	310065
Plastic seat for Neptun	312030
Telescopic legs (1 pc. with clips, rubber feet, bush)	312033
Rubber feet (for telescopic legs)	310017
Bush (for telescopic legs)	310021
Plastic part for backrest	312031

CE

This product is CE marking which means the product complies with EU legislations and enables the free movement of the product in the European market.