

Adjustable Shower Stool



↓As acc. You can get a soft foam cushion. It comes in two widths.



Solid and very stable shower stool, which have been equipped with telescopic legs.

The telescopic legs make it possible to adjust the shower seat from 420 mm to 625 mm in height with 8 different adjustments. This can be done without any tools.

The structure of the telescopic legs and the non-slip rubber feet make sure the shower stool stands very stable.

The shower stool is stackable and is prepared for the mounting of a back rest (easily put on with the help of 4 screws). The back is made of the same strong materials.

You can get the adjustable shower stool in two different widths.





	Technic	al Data:	
	Adjustable s	shower stool	
	Adjustable shower stool		Adjustable shower stool
	Width 410 mm		Width 500 mm
Max. Seat height	625 mm		625 mm
Min. seat height	420 mm		420 mm
Max. Total height	770 mm		770 mm
Min. total height	560 mm		560 mm
Seat width	410 mm		500 mm
Max. Width at floor	595 mm		680 mm
Min. width at floor	560 mm		640 mm
Seat depth	360 mm		360 mm
Max. Depth at floor	520 mm		520 mm
Min. depth at floor	440	mm	440 mm
Weight	4,7	' kg	5,0 kg
Max. Load	150) kg	150 kg
		Rilsan which is a nylon coating with following characteristics: Weatherproof, Chemically proof, Proo against salt mist for over 2000 hours, Physiological harmless, 220 - 250 my	
Slats		Plastic slate 60 mm width and 6 mm thick	
Rubber Feets		Non slip technical grey rubber "no marks at the floor" Hole in the rubber feet for outflow of water	
Telescopic legs		White lacquered aluminium tube fasten with same typ of clips as for elbow crutch	
Shower stool back rest		Max. Width 410 mm, height over the seat 320 mm Mounted with 4 6 mm screws	
Description:		Art. no	
Adjustable shower stool 500 mm		310122	
Adjustable shower stool 410 mm		310120	
Shower stool back rest		310180	
Foam cushion. 36x37 cm (for 41 cm)		310178	
Foam cushion 36x47 cm (for 50 cm)		310179	
Slat middle 410 mm		310151	
Slat edge 410 mm		310156	
Slat middle 500 mm		310150	
Slat edge 500 mm		310155	
Rubber feet		310175	
Screws for shower stool back		310184	
Spacers for shower stool back		310182	

This product is CE marking which means the product complies with EU legislations and enables the free movement of the product in the European market.

07.10.15

