

Neptun Small Shower Stool. Stainless Steel



Small shower stool made for places where you do not have enough room for a standard shower stool.

Made in Stainless Steel (telescopic legs in aluminium).

Very comfortable seat.

Very easy to clean.

Very easy to adjust height (43-58 cm).



⇒Air filled plastic seat which you can click off (for cleaning a.s.o.). The seat is with a nobly surface and with water grooves.

←The frame is designed to make it very easy to clean.



 \Rightarrow The Height is adjusted with clips



NB! You cannot use a backrest with this item.



Technical data: Neptun Small shower Stool. Stainless Steel	
Seat size. (size of the plastic seat)	W41xD36, 5xH4,5 cm
App. Widths at floor	40-42,5 cm
App. Depths at floor	44-46 cm
Weight	3,7 kg
Seat plastic part material	Polythene plastic
Cleaning	Ordinary cleaning detergents at max. 80 degrees.
Stackable	Yes
Max. user load	150 kg
Frame	Stainless Steel.
	Powder lacquered
Telescopic legs	Aluminium. Height adjustable in 8 positions.
Rubber feet	Grey plastic. Hole in bottom for outflow of water.
NB! You cannot use a b	ackrest with this item
Description:	Art. no
Neptun small shower stool	312100
Spare parts:	
Clips for height adjustment	310065
Plastic seat for Neptun	312030
Telescopic legs (1 pc. with clips, rubber feet, bush)	312130
Rubber feet (for telescopic legs)	310017
Bush (for telescopic legs)	312131

CE

This product is CE marking which means the product complies with EU legislations and enables the free movement of the product in the European market.

07.10.15

